



sartorius

Sartorius Liquid Handling Products



turning science **into solutions**

Fapa s.a.s.
Via Pascoletto, 20
24040 Lallio (BG)
Tel. 035.6221219
Fax. 035.4372675
fapa@fapa.bg.it
www.fapa.bg.it



|| The Sartorius liquid handling products are a perfect combination of ergonomics, reliability and design.

■ Table of Contents

□ About Sartorius

- 7 About Us
- 9 Responsible Manufacturing
- 10 Outstanding Product Quality
- 12 Ergonomics, Reliability and Design

□ Pipette Selection Guide

- 15 Electronic or Mechanical Pipette

□ Electronic Pipettes

- 18 Picus
- 26 eLINE®

□ Mechanical Pipettes

- 34 mLINE®
- 40 Proline® Plus
- 46 Proline®

□ Pipette Tips

- 52 Optifit Tips and SafetySpace™ Filter Tips
- 54 SafetySpace™ Filter Tips
- 56 Package Options
- 62 Tip Compatibility Charts

□ Dispensers and Tips

- 68 eLINE® Lite and Pro Dispensers
- 70 Dispenser Tips
- 72 Mechanical Stepper

□ Stands

- 74 Pipette Stands

□ Accessories

- 80 Safe-Cone Filters
- 82 Elbow Pad
- 83 Reagent Vessel
- 83 Cooling Rack

□ Maxi-volume Liquid Handling

- 86 Midi Plus™
- 88 Biofiller
- 90 Proline® Prospenser
- 91 Prospenser
- 92 Biotrate Digital Burettes

□ Pipetting Academy

- 96 Pipetting Academy™
- 98 Pipetting Recommendations

□ Pipette Calibration and Maintenance Services

- 102 Pipette Calibration and Maintenance Services
- 108 Quanta – Pipette Service and Calibration Software
- 112 Pipette Decontamination Procedure
- 113 Autoclaving Instructions
- 114 Troubleshooting Guide




sartorius

Proline Plus

10.00

0.5 - 10 µl



■ About Us


Sartorius Biohit Liquid Handling, a part of Sartorius Group, is a leading, global provider of electronic and mechanical pipettes, disposable pipette tips and related services, for pharmaceutical and chemical industry, clinical laboratories, research institutes and universities.

Known as the forerunner in developing ergonomic and light pipettes for years with the users comfort and health in mind, Sartorius is a trusted supplier to many laboratory professionals who want to reduce their risk of work-related injuries. Strong technical innovations and the use of the newest technologies both in design and manufacturing ensure the highest reliability and quality of Sartorius' liquid handling products. Sartorius offers pipette maintenance, repair and calibration services globally through its certified service centres.

Sartorius' liquid handling business is based in Helsinki, Finland. Its own manufacturing sites are both in Finland and in China and the sales network is global, shared with the Sartorius Group, in 110 countries.

Sartorius is one of the world's leading providers of laboratory and process technologies and equipment covering the segments of Bioprocess Solutions, Lab Products & Services and Industrial Weighing. Founded in 1870, the Gottingen, Germany, based company currently employs more than 5,500 people around the world. Sartorius has its own production facilities in Europe, Asia and America.





|| We are committed to our environmental policy that dictates the rules for our product development, manufacturing and logistics.

■ Responsible Manufacturing



We pay particular attention to the environmental impact of our operations. We aim to develop and manufacture products that will cause the smallest possible environmental load throughout their lifecycles.

Environmentally friendly design and manufacturing

Sartorius complies with ISO 14001 environmental standards. Already at the design stage, we look into ways of reducing the usage of hazardous substances and materials. Our products are also designed to minimise waste during manufacturing and use. For example, the tip refill system can reduce waste by up to 61% compared to using racked tips in laboratories. In production, we have invested in technologies that generate less waste. We also use environmentally friendly carbon-free electricity in our production facility in Kajaani, Finland.

Recyclable materials

The plastic materials used for the products and their packaging are chosen to be suitable for use in waste-to-energy facilities as far as possible. For example the tips and the tip racks made of 100% polypropylene (PP) can be fully recycled as energy waste – PP could also be reused. The cardboard used for packaging can be recycled (reused).

□ In Practice

- Package materials are suitable for recycling either for reuse or as energy waste
- Package sizes are minimised in order to use less material and to render logistics more efficient
- Pipette tips and racks are 100% recyclable as energy waste or for reuse
- Cadmium-free plastics have been used since 1994
- Cadmium-free batteries have been used since 1996
- mLINE and Proline Plus mechanical pipettes are more than 90% recyclable as energy waste
- Electronic pipettes are WEEE/RoHS compliant, which means controlled recycling by national authorities
- Environmentally friendly carbon-free electricity is used in our production facility in Kajaani, Finland
- Minimum possible amount of paper and energy is used in our offices
- Environmentally friendly paper is used for printing

■ Outstanding Product Quality

Sartorius' products are developed and manufactured according to the requirements of the ISO 9001, ISO 13485 and ISO 14001 quality and environmental standards. The tip production also follows standard ISO 14644-1 in order to fulfil ISO class 8 cleanroom conditions. Our accredited pipette calibration follows the ISO 17025 standard. Our pipettes are supplied with individual quality control certificates.

We aim to continuously develop our products and processes in order to meet – and often exceed – the demands of regulatory authorities, environmental bodies, and most importantly, our customers.



Sartorius' products are developed and manufactured according to the requirements of the ISO 9001, ISO 13485 and ISO 14001 quality and environmental standards. Tip production also abides by the ISO 14644-1 standard, in order to fulfil ISO Class 8 Cleanroom conditions. ISO 13485 is a specific standard for medical device quality systems, and supplementing the more generic ISO 9001 standard, which applies to many industries.



The Finnish national accreditation body operates independently as part of the Measurement Technology Centre (MIKES). Accredited pipette calibration laboratories in Finland, Germany, France, UK, China and Japan calibrate pipettes according to precise technical requirements. Our calibration laboratories in Finland, Germany, France, UK, China and Japan have been granted this status by their national accreditation bodies.



Sartorius offers a 2-year warranty for all mechanical and electronic pipettes. The low lifetime cost and environmental friendliness of our products, which have long warranty periods, give a high return on investment.

We follow these manufacturing quality standards

ISO 9001 · ISO 13485 · ISO 14001 · ISO 17025



The ergonomic design label indicates products, which Sartorius has specifically designed to reduce the risk of work-related hand, arm and shoulder disorders, such as RSI (Repetitive Strain Injury).



The Optiload tip loading mechanism developed by Sartorius in mLINE, Proline Plus, eLINE and Picus pipettes allows tips to be loaded with constant force. In turn, this secures optimal tip sealing and minimum tip ejection force.



DNase, RNase and endotoxin free certification is issued for each lot of Sartorius Single Tray and Refill Pack tips, in order to protect the sample from contamination. This certificate can be downloaded from www.biohit.com. Sartorius' tip production is ISO 8 cleanroom classified, in order to secure contamination free manufacturing environment and products.



Most Sartorius liquid handling products are autoclavable. Please see details in product descriptions.

■ Ergonomics, Reliability and Design

Three key factors – ergonomics, reliability and design – form the cornerstone of all of our R&D projects for new products. The newest family member, the Picus electronic pipette, is an excellent example of combining all of these aspects in one product: it is the lightest and smallest electronic pipette on the market. As well as winning a Red Dot design and a Fennia Prize award, this device is highly reliable, accurate and precise. All Sartorius Liquid Handling products are designed in Finland, where our R&D team is constantly seeking solutions that make the pipetting experience even better.

“ Designing products that people work with on a daily basis is always challenging. Many users are interviewed and multiple aspects need to be taken into account, to combine excellent ergonomics and easy usability with today's technology and features. To solve this puzzle and come up with a great product is an exciting, but sometimes tough, journey. However, it is always rewarding in the end.”

– Ville Hintikka, Chief Designer at Sartorius



□ Ergonomics

When designing a pipette, we always consider the shape and function of the human hand. Because we understand the risks of repetitive pipetting, we emphasise ergonomic design in every product we make. Simply put, this means products that you can use in a comfortable posture with minimum muscle power. Our pipettes and dispensers are designed for both right- and left-handed users. Their operating buttons are located sufficiently close together, within ergonomic reach of the thumb.



□ Reliability

For us, reliability has many aspects, the most important of which are the accuracy and precision of the results and the secured purity.

The core of a pipette lies in its **accuracy and precision**. For this reason, we have used the newest technologies together with in-house innovations, to achieve even more reliable pipetting results. Our electronic brake, piston control system and plate tracker for electronic pipettes are our newest innovations. They increase the devices' accuracy, precision and reliability. Another important factor in achieving reliable results is the optimal tip fit, which we can guarantee by designing and producing the tips ourselves, to perfectly match our pipettes.



Because **purity** is a key concern in many laboratories, we offer special Safe-Cone Filters for our pipettes, to keep them clean. We strive to produce as many of our products as possible, pipettes and tips, to be autoclavable. Our pipette tips are manufactured in ISO Class 8 Cleanroom conditions and we test every certified tip lot, for DNase, RNase and endotoxins, in an external laboratory. We also offer an innovative SafetySpace™ filter tip range for safer and contamination-free pipetting.

□ Design



reddot design award
product design

We aim to provide products with a timeless and light, yet practical, design, suitable for laboratory settings and pleasing to the eye of the user. The newest member of our product family, the Picus electronic pipette, won the Red Dot design award and the Fennia Prize Honorary Mention in 2012. Our other electronic pipette family, eLINE, was given an honorary mention in the Pro Finnish Design competition at the time of its launch.





0.2 - 10 µl



■ Pipette Selection Guide

□ Electronic or Mechanical Pipette

Are you looking for a tool for sterile work, or for a pipette you could easily calibrate yourself? Or are you looking for a particularly comfortable solution, something really light and ergonomic? By consulting the tables below, you can choose the tool that best suits you and your work.

Electronic or Mechanical Pipette

Features	Electronic Pipettes	Mechanical Pipettes
Highest ergonomics	•	
Fastest pipetting	•	
User-independent results	•	
Multiple pipetting modes	•	
Fully autoclavable		•
Calibration by user	• ¹	•

Electronic Pipettes

Features	Picus	eLINE
Most ergonomic	•	
Weight ²	100 g	170 g
Multiple pipetting modes ³	8+5	7+3
Microwell plate tracker	•	
Electronic tip ejection	•	•
Calibration by user	•	
Information on service & calibration intervals	•	
Hot key for stored programs	•	
Memory places	10	6
Safe-Cone Filters	•	•
Autoclavable lower parts	•	•
Optiload in multichannels	•	•
Colour coding on pipette	•	•
Warranty for 2 years	•	•

Mechanical Pipettes

Features	mLINE	Proline Plus	Proline
Most ergonomic	•		
Pipetting forces ⁴	12 N	15 N	20 N
Fully autoclavable	•	•	
Optiload tip loading mechanism	all models	multichannels only	
Weight ⁴	77 g	82 g	84 g
Safe-Cone Filters	•	•	•
Filter ejector	•		
Volume locking	•	click stops	click stops
Thermal insulation	•	•	
Colour coding on pipette	•	•	
Colour coding caps	•		•
ID tags	•		
Warranty for 3 years	•	•	•

1) Picus only

2) 300 µl 1-channel models

3) Include mode additions

4) 1000 µl 1-channel models



MENU EDIT ADV

Pipetting

300.0 μ l

7

Picus
300

■ Electronic Pipettes

Table of Contents

18 Picus

26 eLINE®

■ Picus Electronic Pipettes

The most sophisticated and ergonomic pipette ever!

Picus, the winner of the Red Dot design award and the Fennia Prize Honorary Mention in 2012, is Sartorius' newest, and most ergonomic, pipette in the Biohit*family*. Designed to revolutionise pipetting, it is an extremely small and lightest electronic pipette on the market. This enables it to ease the user's workload and provides protection from repetitive strain injury (RSI). Its new generation technology, electronic brake and piston control system guarantee accurate and precise pipetting results. By guiding pipetting steps, the unique plate tracker increases reliability in microwell plate work. Picus is available in single-channel models, covering a volume range of 0.2-10 000 µl and in multichannel models from 0.2 µl to 1200 µl.

- New ergonomic design protects you from RSI (repetitive strain injury) and eases your workload in long pipetting series
- Our lightest and smallest electronic pipette ensures an optimal working posture and offers comfortable pipetting
- New generation technology – electronic brake and piston control system guarantee accurate and precise pipetting results
- Intuitive user interface – learn main functions easily without a manual



reddot design award
product design



fenniaprize 12
good design grows global
HONORARY MENTION



Enjoy Pipetting

- New ergonomic design protects you from RSI (repetitive strain injury) and eases your workload in long pipetting series
- Our lightest and smallest electronic pipette ensures an optimal working posture and comfortable pipetting
- Effortless, electronic tip ejection

Achieve Reliable Results

- New generation technology – the electronic brake and piston control system guarantee accurate and precise pipetting results
- A unique plate tracker increases reliability in microwell plate work, by guiding pipetting steps
- Fully electronic operation guarantees user-independent results

Learn to Use it in a Minute!

- Learn the main functions easily without a manual
- Fastest volume setting and easy menu navigation with unique adjustment wheel
- Hot key – fast access to your favourite programs, easing routine work

□ Features and Benefits

Hot key for fast access to up to 10 saved programs

Soft keys for programming

Unique electronic tip ejection

Operating button

Adjustment wheel

- Offers extremely fast volume setting and easy menu navigation
- Enables ergonomic one-handed operation
- Controls accurately manual pipetting and titrating speed with just a light touch of the thumb


Intuitive user-interface

Ergonomics

- Extremely small and light electronic pipette (only 100 g)
- Comfortable handle design and finger hook
- Unique electronic tip ejection
- Together with our highly regarded Pipetting Academy training programme, Picus will help you to improve your health and safety in the lab



Autoclavable lower parts (except for 1200 µl multi-channel lower parts)



Optiload for perfect tip sealing

- Allows tip loading with an equal constant force to every channel
- Enables perfect tip sealing to every individual tip cone



Safe-Cone Filters in all models > 10 µl to prevent contamination






New generation technology

- Enhanced DC-motor concept takes Picus a step beyond existing technological solutions, providing outstanding performance
- Electronic brake stops piston movement rapidly and precisely, ensuring high precision, especially in serial dispensing
- Optical sensor controls and monitors piston movement in real time, ensuring unbeatable accuracy and reliability



Microwell plate tracker

- Unique built-in tracker for 96 and 384 well plates guides the user to pipette into the correct wells
- Improves work efficiency and reliability of results
- Tracker is an optional feature used in conjunction with Pipetting, Reverse Pipetting and Multi Dispensing modes



Technical Features

- Calibration adjustment in 1, 2 or 3 points
- Information for service and calibration intervals
- Autoclavable lower parts (excl. 1200 µl models)
- Charging in charging stand or with micro USB cable
- Possible to continue working with USB charging on
- Li-Polymer battery enables charging time of approx. 1 hour



Pipetting Modes

**Main Mode-
Available in All Picus Models**

Additional Functions - to be Used in Conjunction with the Main Mode

		Tracker	Mixing	Counter	Excess Volume Adjustment	Auto Dispensing
Pipetting	✓	✓	✓	✓		
Reverse Pipetting	✓	✓		✓	✓	
Manual Pipetting	✓					
Multi Dispensing	✓	✓			✓	✓
Diluting	✓		✓			
Sequential Dispensing	✓				✓	
Multi Aspiration	✓					
Titrate	✓					

Ordering Information


Picus

Cat. No.	Channels	Volume Range (µl)	Increment (µl)	Test Volume (µl)	Inacc. (%)	Impr. (%)	Safe-Cone Standard	Filters Plus	Optifit Tips (µl)	SafetySpace™ Filter Tips (µl)
735021	1	● 0.2-10	0.01	10	0.90	0.40	-	-	10	10
				5	1.00	0.70				
				1	2.50	1.50				
735041	1	● 5-120	0.10	120	0.40	0.15	721008	721018	200 350	120
				60	0.60	0.20				
				12	2.00	1.00				
735061	1	● 10-300	0.20	300	0.40	0.15	721007	721017	350	300
				150	0.60	0.20				
				30	1.50	0.80				
735081	1	● 50-1000	1.00	1000	0.40	0.15	721006	721016	1000	1000
				500	0.60	0.20				
				100	1.50	0.50				
735101	1	● 100-5000	5.00	5000	0.50	0.15	721005	721015	5000	-
				2500	0.80	0.20				
				500	1.00	0.40				
735111	1	● 500-10000	10.00	10000	0.60	0.20	721005	721015	10000	-
				5000	1.20	0.30				
				1000	3.00	0.60				
735321	8	● 0.2-10	0.01	10	0.90	0.50	-	-	10	10
				5	1.50	0.80				
				1	4.00	3.00				
735341	8	● 5-120	0.10	120	0.50	0.20	721008	721018	200 350	120
				60	0.70	0.30				
				12	2.00	1.50				
735361	8	● 10-300	0.20	300	0.50	0.20	721007	721017	350	300
				150	0.70	0.30				
				30	2.00	1.00				
735391	8	● 50-1200	1.00	1200	0.50	0.20	721006	721016	1200	1200
				600	1.00	0.30				
				120	2.50	1.00				
735421	12	● 0.2-10	0.01	10	0.90	0.50	-	-	10	10
				5	1.50	0.80				
				1	4.00	3.00				
735441	12	● 5-120	0.10	120	0.50	0.20	721008	721018	200 350	120
				60	0.70	0.30				
				12	2.00	1.50				
735461	12	● 10-300	0.20	300	0.50	0.20	721007	721017	350	300
				150	0.70	0.30				
				30	2.00	1.00				
735491	12	● 50-1200	1.00	1200	0.50	0.20	721006	721016	1200	1200
				600	1.00	0.30				
				120	2.50	1.00				

All pipettes include universal AC-adaptor (EU, UK, US | JPN, AUS and CHN plugs)



Charging Stand



Charging Carousel



Linear Stand

□ Picus Charging Options

- 1-Place Charging Stand
- 4-Place Charging Carousel
- USB cable

Picus electronic pipettes can be charged with a 1-place charging stand, a 4-place charging carousel or a direct charging USB cable. The variety of charging options has been designed for optimal utility.

It is possible to continue pipetting while the Picus is being charged, using a USB cable. The compact design of the 4-pipette charging carousel is ideal for saving bench space in the laboratory. The 4-place rotating head provides easy access to the desired unit.

□ Ordering Information

Pipette Stands

Cat. No.	Item
730981	Charging Stand for one pipette
730991	Charging Carousel for 4 pipettes
725620	Linear Stand (non-charging)

Charging stands include universal AC-adaptor (EU, UK, US | JPN and CHN plugs)



■ eLINE® Electronic Pipettes


Proven Performance and User Comfort


The electronic Sartorius eLINE pipettes in the Biohit *family* have been providing high precision and comfort for years. This award-winning (ProFinnish Design 2001) device is the first electronic pipette designed to take care of its user. The comprehensive range of pipetting modes in eLINE reduces the need for several work stages and enables liquid dispensing twice as fast as a mechanical pipette. Its unique DC-motor concept, with build-in error control, improves pipetting precision and allows more reliable results. eLINE is available as single-channel models, covering a volume range of 0.1 to 5000 µl, and as multichannel models from 0.2 to 1200 µl.

- Proven accuracy and precision with DC-motor concept
- Fully electronic operations guarantee user-independent results and prevent RSI
- Effortless, electronic tip ejection
- Twice as fast as a mechanical pipette
- One eLINE covers the volume range of two mechanical pipettes
 - fewer pipettes needed
- Comprehensive range of pipetting modes (8)
- Several charging options
- Liquid dispensing from the air, without touching the vessel wall, with eLINE 0.1 – 5 µl model



□ Features and Benefits





The increased liquid dispensing speed of the eLINE 0.1–5 µl pipette with Super Pipetting feature enables liquid dispensing from the air, emptying the tip completely, without hanging last drop.

Contamination Free Pipetting with Super Pipetting (SP) Feature

It is now possible to dispense the very smallest volumes from the air – without needing to touch the receiving vessel wall or for dispensing under the liquid surface. The increased liquid dispensing force of the SP feature empties the tip completely, with no drops left hanging on to the tip. Using the SP feature, the user can also prevent droplets from clinging onto the electrically charged outer surface of a microplate well, since there is no need to touch the well wall. In both cases, the accuracy of the dispensed volume increases.

- Reduced risk of contamination
- Improved accuracy for small volumes
- Available in eLINE 0.1–5 µl model

Proven Accuracy and Precision

A unique DC-motor concept for 10 times more precise dispensing compared to a stepper-motor-driven system.

- Piston movement is halted using an extremely rapid solenoid brake
- Major advantage in multi-dispensing where the liquid column has to be cut sharply to achieve even doses
- Error control with optical sensor
- Controls and monitors the piston movement in real time
- Notifies the user if the piston is in an incorrect position

Pipetting Modes

Main Mode

Additional Functions

		Mixing	Counter	Auto Dispensing (Timed)
Pipetting	✓	✓	✓	
Reverse Pipetting	✓			
Manual Pipetting ¹	✓			
Multi Dispensing	✓			✓
Diluting	✓	✓		
Sequential Dispensing ²	✓			
Multi-Aspirating	✓			
Super Pipetting	✓	Only available in eLINE 0.1–5 µl		

1) Not available in eLINE multichannel pipettes

2) Not available in eLINE 0.1–5 µl

Ordering Information

eLINE®


Cat. No.	Channels	Volume Range (µl)	Increment (µl)	Test Volume (µl)	Inacc. (%)	Impr. (%)	Safe-Cone Standard	Filters Plus	Optifit Tips (µl)	SafetySpace™ Filter Tips (µl)
73001x	1	● 0.1-5	0.05	5	1.20	0.70	-	-	10	10
				2.5	1.70	1.20				
				0.5	9.50	8.00				
73002x	1	● 0.2-10	0.05	10	0.90	0.40	-	-	10	10
				5	1.00	0.70				
				1	3.50	2.30				
73004x	1	● 5-120	0.50	120	0.40	0.15	721008	721018	200 350	120
				60	0.60	0.20				
				12	2.00	1.00				
73006x	1	● 10-300	1.00	300	0.40	0.15	721007	721017	350	300
				150	0.60	0.20				
				30	1.50	0.80				
73008x	1	● 50-1000	5.00	1000	0.40	0.15	721006	721016	1000	1000
				500	0.60	0.20				
				100	1.50	0.50				
73010x*	1	● 100-5000	10.0	5000	0.50	0.15	721006	721016	5000	-
				2500	0.80	0.20				
				500	1.00	0.40				
73032x	8	● 0.2-10	0.05	10	0.90	0.50	-	-	10	10
				5	1.50	0.80				
				1	4.00	3.00				
73034x	8	● 5-120	0.50	120	0.80	0.20	721008	721018	200 350	120
				60	0.70	0.30				
				12	3.00	1.50				
73036x	8	● 10-300	1.00	300	0.50	0.20	721007	721017	350	300
				150	0.70	0.30				
				30	2.00	1.00				
73039x	8	● 50-1200	5.00	1200	0.50	0.20	721006	721016	1200	1200
				600	1.00	0.30				
				120	2.50	1.00				
73042x	12	● 0.2-10	0.05	10	0.90	0.50	-	-	10	10
				5	1.50	0.80				
				1	4.00	3.00				
73044x	12	● 5-120	0.50	120	0.80	0.20	721008	721018	200 350	300
				60	0.70	0.30				
				12	3.00	1.50				
73046x	12	● 10-300	1.00	300	0.50	0.20	721007	721017	350	300
				150	0.70	0.30				
				30	2.00	1.00				
73049x	12	● 50-1200	5.00	1200	0.50	0.20	721006	721016	1200	1200
				600	1.00	0.30				
				120	2.50	1.00				

x = 0 without AC-adaptor, x = 1 with AC-adaptor (EU, UK, US | JPN and CHN plugs)

* NOTE: Min. volume in P-mode is 500 µl. 100 µl is possible in d-mode



Charging Stand



Charging Carousel



Linear Stand

□ eLINE® Charging Options

- 1-Place Charging Stand
- 4-Place Charging Carousel
- AC-adaptor

The eLINE electronic pipettes can be charged with a 1-place charging stand, a 4-place charging carousel or a direct charging AC-adaptor. The variety of charging options has been designed for optimal utility. It is possible to

continue pipetting while the eLINE is being charged through the AC-adaptor. The compact design of the 4-pipette charging carousel is ideal for saving bench space in a laboratory. The 4-place rotating head provides easy access to the desired unit.

□ Ordering Information

Pipette Stands and Accessories

Cat. No.	Item
730981	Charging Stand for one pipette
730991	Charging Carousel for 4 pipettes
725620	Linear Stand (non-charging)
731009	eLINE replacement battery

Charging stands include universal AC-adaptor (EU, UK, US | JPN and CHN plugs)





sartorius

sartorius

mLINE
200
20 - 200 µl

mLINE
100.0
10 - 100 µl

100 - 1000 µl

■ Mechanical Pipettes

Table of Contents

34 mLINE®

40 Proline® Plus

46 Proline®

■ mLINE® Mechanical Pipettes

Effortless Accuracy

Sartorius' most advanced mechanical pipette in the Biohit*family* – the mLINE – offers excellent ergonomics, performance and safety in manual pipetting. It is particularly designed for repetitive and long-lasting pipetting in order to prevent work related hand, arm and shoulder disorders or Repetitive Strain Injury (RSI). mLINE is recommended by health and safety officers around the globe. It covers the full volume range of 0.1 µl to 10 ml and is available in single- and multichannel models.

- Excellent ergonomics to protect from injuries, with the exceptionally light pipetting and tip ejection forces
- High accuracy and precision, also in repetitive, long-lasting pipetting
- Minimised risk of contamination with Safe-Cone Filters and full autoclavability
- Increased safety with volume lock which prevents accidental volume changes while pipetting
- Thermal insulation of internal components to improve accuracy



Extremely light pipetting forces & ergonomic design

- Prevents repetitive strain injury (RSI)
- Better pipetting results in long pipetting series

Ergonomic finger support

- No need to squeeze the pipette during pipetting, reduces the risk of RSI



□ Features and Benefits

Smooth and light volume setting
– Prevents hand injuries, can be adjusted by finger tips

Thermal insulation of internal components to improve accuracy

Volume colour-coding
– Eases correct Sartorius tip selection, through matching colour code

Unique thumb-activated light filter-ejection mechanism

Materials have high chemical and UV-resistance
– Highly durable and long lasting pipette

Large numbers used on volume display
– Quick and easy-to-read volume minimises risk of volume-setting errors

Volume lock
– Prevents accidental volume changes while pipetting

Easy to use "do-it-yourself" calibration
– Quick re-calibration possible e.g. for different liquids

Only three parts to disassemble
– The easiest pipette to clean and maintain, improves purity

Fully autoclavable without disassembly (excl. 1200 µl multi-channel pipettes)
– Easy to steam sterilise, improves purity



Fully autoclavable

Optiload for perfect tip sealing
– Allows tip loading with an equal constant force onto every channel
– Enables perfect tip sealing onto every individual tip cone
– Makes tip ejection light and easy



Safe-Cone Filters available for models > 10 µl
– Contamination-free pipetting





Superior Ergonomics

mLINE pipettes have a patented spring mechanism (plunger technology) resulting in exceptionally light pipetting forces. The tension spring mechanism makes the pipetting force almost 4 times lighter than conventional compression spring operated pipettes. Light pipetting force improves pipetting precision in long pipetting series. Thanks to the device's unique spring mechanism, the starting force is always constant, regardless of the set volume, which improves pipetting results especially for small volumes.



Perfect Tip Sealing with Optiload Tip Loading Mechanism

All mLINE pipettes are equipped with spring loaded tip cones – the Optiload mechanism secures even tip sealing and simultaneously allows tips to be loaded and ejected with minimum force. This is a particular advantage with multichannel models, where tip loading and ejecting generally requires more force. Using less force reduces the risk of hand injuries.



Prevent Contamination with Safe-Cone Filters

Replaceable Safe-Cone Filters act as a barrier, preventing aerosols or fluids from reaching the internal components of the pipette. Safe-Cone Filters are available for all mLINE models greater than 10 µl. These filters must be replaced regularly, and in every case of over-aspiration.

Unique Safe-Cone Filter Ejection Mechanism

Safe-Cone filters can easily and safely be ejected by removing the colour cap and pressing all the way down on the operating button.

Ordering Information

mLINE®

Cat. No.	Channels	Volume Range (µl)	Increment (µl)	Test Volume (µl)	Inacc. (%)	Impr. (%)	Safe-Cone Standard	Filters Plus	Optifit Tips (µl)	SafetySpace™ Filter Tips (µl)
725010	1	● 0.1-3	0.002	3	1.30	0.80	-	-	10	10
				1.5	2.40	1.60				
				0.3	10.00	6.00				
725020	1	● 0.5-10	0.01	10	1.00	0.60	-	-	10	10
				5	1.50	1.00				
				1	2.50	1.50				
725030	1	● 2-20	0.02	20	0.90	0.40	721014	-	200	20
				10	1.20	1.00				
				2	3.00	2.00				
725050	1	● 10-100	0.10	100	0.80	0.15	721008	721018	200	120
				50	1.00	0.40		350		
				10	2.00	1.00				
725060	1	● 20-200	0.20	200	0.60	0.15	721007	721017	200	200
				100	0.80	0.30		350	300	
				20	2.00	0.80				
725070	1	● 100-1000	1.00	1000	0.60	0.20	721006	721016	1000	1000
				500	0.60	0.20				
				100	1.00	0.40				
725080	1	● 500-5000	10.0	5000	0.50	0.20	721005	721015	5000	-
				2500	0.60	0.30				
				500	2.00	0.60				
725090	1	● 1-10 ml	20.0	10000	0.60	0.20	721005	721015	10 000	-
				5000	1.20	0.30				
				1000	3.00	0.60				
725120	8	● 0.5-10	0.01	10	1.50	1.00	-	-	10	10
				5	2.50	2.50				
				1	4.00	4.00				
725130	8	● 5-100	0.10	100	0.70	0.25	721008	721018	200	120
				50	1.00	0.70		350		
				10	3.00	1.50				
725140	8	● 30-300	0.20	300	0.60	0.25	721007	721017	350	300
				150	1.00	0.50				
				30	2.00	1.00				
725220	12	● 0.5-10	0.01	10	1.50	1.00	-	-	10	10
				5	2.50	2.50				
				1	4.00	4.00				
725230	12	● 5-100	0.10	100	0.70	0.25	721008	721018	200	120
				50	1.00	0.70		350		
				10	3.00	1.50				
725240	12	● 30-300	0.20	300	0.60	0.25	721007	721017	350	300
				150	1.00	0.50				
				30	2.00	1.00				


□ Ordering Information

Pipette Stands and Accessories

Cat. No.	Item
725600	Carousel Stand for 6 pipettes
725610	Pipette Holder for one pipette
725620	Linear Stand
726203	Calibration tool Tube opener
726001	Colour coding caps, 5 pcs



Carousel Stand



Pipette Holder



Linear Stand



Calibration tool





□ mLINE® Starter Kits and mLINE® PCR Starter Kit


mLINE Starter Kit offers an opportunity to test and get started with mLINE pipettes. You can choose between two Starter Kits with three single-channel pipettes of different volumes. Both kits also include colour coded tip racks for every pipette and many useful accessories such as a pipette holder,

colour coding caps and a calibration tool, which also acts as a tube opener. mLINE PCR Starter Kit includes all you need to make your PCR work fast and reliable. A cooling rack keeps valuable samples stable while working, thereby increasing reliability.



□ Ordering Information

Cat. No.	Item
725651	mLINE® Starter Kit 1 mLINE 0.5 – 10 µl mLINE 10 – 100 µl mLINE 100 – 1000 µl Optifit Tip 0.1 – 10 µl, Single Tray Optifit Tip 0.5 – 200 µl, Single Tray Optifit Tip 10 – 1000 µl, Single Tray Pipette Holder × 3 Colour Coding Caps Calibration Tool Tube Opener × 3 Pen Literature
725652	mLINE® Starter Kit 2 mLINE 2 – 20 µl mLINE 20 – 200 µl mLINE 100 – 1000 µl Optifit Tip 0.5 – 200 µl, Single Tray Optifit Tip 10 – 1000 µl, Single Tray Pipette Holder × 3 Colour Coding Caps Calibration Tool Tube Opener × 3 Pen Literature
725660	mLINE® PCR Starter Kit 3 mLINE 0.5 – 10 µl mLINE 10 – 100 µl mLINE 100 – 1000 µl SafetySpace Filter Tip 0.1 – 10 µl, Single Tray, sterile SafetySpace Filter Tip 2 – 120 µl, Single Tray, sterile SafetySpace Filter Tip 50 – 1000 µl, Single Tray, sterile Pipette Holder × 3 Colour Coding Caps Calibration Tool Tube Opener × 3 Cryo Pen Cooling Rack for Micro Tubes



■ Proline® Plus Mechanical Pipettes

Dependable Durability

Sartorius' mechanical pipette Proline Plus in the Biohit*family* is designed to offer comfort and quality in every day manual pipetting. While it shares many of the excellent features of mLINE, it has a personal design and a robust feel for heavier use. In addition, it has the widest pipette range, including fixed volume pipettes, for when volumes need to be ready-set to avoid errors. Proline Plus pipettes are an excellent choice for both experienced laboratory professionals and students.

- The widest range of pipettes for various users and applications
- Ergonomic design with light pipetting forces, comfortable handle and finger support for reduced risk of strain injuries
- Minimised risk of contamination with Safe-Cone Filters and full autoclavability
- Highly durable with strengthened structure also for heavier use

Widest range of pipettes, both adjustable and fixed volume, available

- Adjustable volume pipettes for professionals
- Fixed volume pipettes for specific applications to avoid volume errors



Optiload tip loading mechanism in multichannel models

- Allows tip loading with an equal constant force onto every channel
- Enables perfect tip sealing onto every individual tip cone
- Makes tip ejection light and easy



□ Features and Benefits

Volume setting with click stop
– Easy to use volume setting

Ergonomic finger support
– No need to squeeze the pipette during pipetting, reduces the risk of RSI

Light pipetting forces & ergonomic design
– Prevents repetitive strain injury (RSI)

Volume colour-coding
– Eases correct tip volume selection through matching colour code

Thermal insulation of internal components
– Improves accuracy

Large numbers in volume display
– Quick and easy-to-read volume minimises risk of volume setting errors

Easy to use "do-it-yourself" calibration
– Quick re-calibration possible e.g. for different liquids

Materials have high chemical and UV-resistance
– Highly durable and long lasting pipette

Only three parts to disassemble
– The easiest pipette to clean and maintain, improves purity

Fully autoclavable without disassembling
– Easy to steam sterilise, improves purity

Safe-Cone Filters available for models > 10 µl
– Contamination free pipetting

Fully autoclavable
121 °C

Ergonomic Design
Sartorius



Widest Range of Pipettes

The Proline Plus family includes the widest range of pipettes, from 3 µl | 5µl to 10 ml, as both adjustable and fixed volume devices. The fixed volume pipettes are particularly useful in student laboratories or in any laboratory where the application requires a constant volume, ready-set, to avoid errors.

Ergonomic Design

A comfortable handle and ergonomic finger support enable effortless pipetting. There is no need to squeeze the handle during pipetting, thereby reducing the risk of RSI.



Perfect Tip Sealing with Optiload Tip Loading Mechanism in Multichannel Pipettes

All Proline Plus multichannel pipettes are equipped with spring loaded tip cones – the Optiload mechanism secures even tip sealing and simultaneously allows tips to be loaded and ejected with minimum force. This is a particular advantage in multichannel models, where tip loading and ejecting generally require more force. Using less force reduces the risk of hand injuries.



Prevent Contamination with Safe-Cone Filters

Replaceable Safe-Cone Filters act as a barrier, preventing aerosols or fluids from reaching the internal components of the pipette. Safe-Cone Filters are available for all Proline Plus models greater than 10 µl. These filters must be replaced regularly, and in every case of over-aspiration.

Maintenance and Calibration Made Easy

No opening tools are needed when cleaning and maintaining Proline Plus pipettes, and only three parts need to be cleaned. These pipettes can be easily calibrated using the tool provided with them.

Ordering Information

Proline® Plus

Cat. No.	Channels	Volume Range (µl)	Increment (µl)	Test Volume (µl)	Inacc. (%)	Impr. (%)	Safe-Cone Standard	Filters Plus	Optifit Tips (µl)	SafetySpace™ Filter Tips (µl)
728010	1	● 0.1–3	0.002	3	1.30	0.80	-	-	10	10
				1.5	2.40	1.60				
				0.3	10.00	6.00				
728020	1	● 0.5–10	0.01	10	1.00	0.60	-	-	10	10
				5	1.50	1.00				
				1	2.50	1.50				
728030	1	● 2–20	0.02	20	0.90	0.40	721014	-	200	20
				10	1.20	1.00				
				2	3.00	2.00				
728040	1	● 5–50	0.10	50	1.00	0.30	721008	721018	200	120
				25	1.40	0.50				
				5	3.00	1.50				
728050	1	● 10–100	0.10	100	0.80	0.15	721008	721018	200	120
				50	1.00	0.40			350	
				10	2.00	1.00				
728060	1	● 20–200	0.20	200	0.60	0.15	721007	721017	200	200
				100	0.80	0.30			350	300
				20	2.00	0.80				
728070	1	● 100–1000	1.00	1000	0.60	0.20	721006	721016	1000	1000
				500	0.60	0.20				
				100	1.00	0.40				
728080	1	● 500–5000	10.0	5000	0.50	0.20	721005	721015	5000	-
				2500	0.60	0.30				
				500	2.00	0.60				
728090	1	● 1–10 ml	20.0	10000	0.60	0.20	721005	721015	10 000	-
				5000	1.20	0.30				
				1000	3.00	0.60				
728120	8	● 0.5–10	0.01	10	1.50	1.00	-	-	10	10
				5	2.50	2.50				
				1	4.00	4.00				
728130	8	● 10–100	0.10	100	0.70	0.25	721008	721018	200	120
				50	1.00	0.70			350	
				10	3.00	1.50				
728140	8	● 30–300	0.20	300	0.60	0.25	721007	721017	350	300
				150	1.00	0.50				
				30	2.00	1.00				
728220	12	● 0.5–10	0.01	10	1.50	1.00	-	-	10	10
				5	2.50	2.50				
				1	4.00	4.00				
728230	12	● 10–100	0.10	100	0.70	0.25	721008	721018	200	120
				50	1.00	0.70			350	
				10	3.00	1.50				
728240	12	● 30–300	0.20	300	0.60	0.25	721007	721017	350	300
				150	1.00	0.50				
				30	2.00	1.00				

Ordering Information

Proline® Plus FIXED Volume, Single-Channel

Cat. No.	Channels	Volume (µl)	Increment (µl)	Test Volume (µl)	Inacc. (%)	Impr. (%)	Safe-Cone Standard	Filters Plus	Optifit Tips (µl)	SafetySpace™ Filter Tips (µl)
728515	1	5	–	5	1.30	1.20	–	–	10	10
728520	1	10	–	10	0.80	0.80	–	–	10	10
728530	1	20	–	20	0.60	0.50	721014	–	200	20
728535	1	25	–	25	0.50	0.30	721008	721018	200	120
728545	1	50	–	50	0.50	0.30	721008	721018	200	120
728550	1	100	–	100	0.50	0.30	721008	721018	200 350	120
728560	1	200	–	200	0.40	0.20	721007	721017	200 350	200 300
728565	1	250	–	250	0.40	0.20	721006	721016	1000	500 1000
728567	1	500	–	500	0.30	0.20	721006	721016	1000	500 1000
728570	1	1000	–	1000	0.30	0.20	721006	721016	1000	1000
728575	1	2000	–	2000	0.30	0.15	721005	721015	5000	–
728580	1	5000	–	5000	0.30	0.15	721005	721015	5000	–
728590	1	10 ml	–	10000	0.60	0.20	721005	721015	10 000	–




Carousel Stand

Ordering Information

Pipette Stands and Accessories

Cat. No.	Item
725600	Carousel Stand for 6 pipettes
725610	Pipette Holder for one pipette
725620	Linear Stand
726203	Calibration Tool Tube opener



Pipette Holder



Linear Stand



Calibration tool



□ Proline® Plus Starter Kits

Proline Plus Starter Kit offers an opportunity to test and get started with Proline Plus. You can choose between four Starter Kits, with two or three single-channel Proline Plus pipettes.

All kits include a range of useful accessories, such as a pipette holder and a calibration tool that also acts as a tube opener, as well as the required instructions.



□ Ordering Information

Cat. No.	Item
728650	Proline® Plus Starter Kit 1 Proline Plus 0.1 – 3 µl Proline Plus 0.5 – 10 µl Optifit Tip 0.1 – 10 µl, Single Tray × 2 Pipette Holder × 2 Calibration Tool Tube Opener × 2 Literature
728651	Proline® Plus Starter Kit 2 Proline Plus 0.5 – 10 µl Proline Plus 10 – 100 µl Proline Plus 100 – 1000 µl Optifit Tip 0.1 – 10 µl, Single Tray Optifit Tip 0.5 – 200 µl, Single Tray Optifit Tip 10 – 1000 µl, Single Tray Pipette Holder × 3 Calibration Tool Tube Opener × 3 Literature
728652	Proline® Plus Starter Kit 3 Proline Plus 2 – 20 µl Proline Plus 20 – 200 µl Proline Plus 100 – 1000 µl Optifit Tip 0.5 – 200 µl, Single Tray × 2 Optifit Tip 10 – 1000 µl, Single Tray Pipette Holder × 3 Calibration Tool Tube Opener × 3 Literature
728653	Proline® Plus Starter Kit 4 Proline Plus 500 – 5000 µl Proline Plus 1000 – 10000 µl Optifit Tip 100 – 5000 µl, Bulk Optifit Tip 1 – 10 ml, Bulk Pipette Holder × 2 Calibration Tool Tube Opener × 2 Literature



Proline® Mechanical Pipettes

Affordable Reliability

The fact that the first mechanical pipette in the Biohit*family*, Proline, is still in use in many laboratories, testifies to its timeless, practical design and reliability. Being the most affordable pipette in Sartorius' range of mechanical pipettes, it is ideal for universities and colleges, or any laboratory seeking a cost-efficient liquid handling tool. Due to its relatively light weight, high accuracy and precision, it is also used by many professionals.

Features and Benefits

Wide range of pipettes, both adjustable and fixed volume

Ergonomic finger support
– No need to squeeze the pipette during pipetting, reduces the risk of RSI


Clear volume display
– Easy to read volume minimises the risk of volume setting errors

"Do-it-yourself" calibration
– Re-calibration possible e.g. for different liquids

Tip cone material highly chemically resistant
– Highly durable tip cone

Safe-Cone Filters available for models > 10 µl
– Contamination-free pipetting

Volume setting with click stop
– Easy-to-use volume setting



Ordering Information

Proline®

Cat. No.	Channels	Volume Range (µl)	Increment (µl)	Test Volume (µl)	Inacc. (%)	Impr. (%)	Safe-Cone Standard	Filters Plus	Optifit Tips (µl)	SafetySpace™ Filter Tips (µl)
720010	1	● 0.1–2.5	0.05	2.5 1.25 0.25	2.50 3.00 12.00	2.00 3.00 6.00	–	–	10	10
720015	1	● 0.5–10	0.10	10 5 1	1.00 1.50 2.50	0.80 1.50 1.50	–	–	10 200	10 20
720080	1	● 2–20	0.50	20 10 2	0.90 1.20 3.00	0.40 1.00 2.00	721008	721018	200	20 120
720025	1	● 5–50	0.50	50 25 5	0.60 0.90 2.00	0.30 0.60 2.00	721008	721018	200 350	120
720050	1	● 10–100	1.00	100 50 10	0.80 1.00 3.00	0.20 0.40 1.00	721007	721017	200 350	120 200 300
720070	1	● 20–200	1.00	200 100 20	0.60 0.80 2.50	0.20 0.30 0.80	721007	721017	200 350	200 300
720060	1	● 100–1000	5.00	1000 500 100	0.60 0.70 2.00	0.20 0.25 0.70	721006	721016	1000	1000
720110	1	● 1000–5000	50.0	5000 2500 1000	0.50 0.60 0.70	0.20 0.30 0.30	721005	721015	5000	–
720210	8	● 0.5–10	0.10	10 5 1	1.50 2.50 4.00	1.50 2.50 4.00	–	–	10	10
720220	8	● 5–50	0.50	50 25 5	1.00 1.50 3.00	0.50 1.00 2.00	721014	–	200 350	120 200 300
720240	8	● 50–300	5.00	300 150 50	0.70 1.00 1.50	0.25 0.50 0.80	721014	–	350	300
720310	12	● 0.5–10	0.10	10 5 1	1.50 2.50 4.00	1.50 2.50 4.00	–	–	10	10
720320	12	● 5–50	0.50	50 25 5	1.00 1.50 3.00	0.50 1.00 2.00	721014	–	200 350	120 200 300
720340	12	● 50–300	5.00	300 150 50	0.70 1.00 1.50	0.25 0.50 0.80	721014	–	350	300

□ Ordering Information

Proline® FIXED Volume

Cat. No.	Channels	Volume (µl)	Increment (µl)	Test Volume (µl)	Inacc. (%)	Impr. (%)	Safe-Cone Standard	Filters Plus	Optifit Tips (µl)	SafetySpace™ Filter Tips (µl)
722001	1	● 5	-	5	1.30	1.20	-	-	10 200	10 20
722004	1	● 10	-	10	0.80	0.80	-	-	10 200	10 20
722010	1	● 20	-	20	0.60	0.50	721008	721018	200 350	20 120
722015	1	● 25	-	25	0.50	0.30	721008	721018	200 350	120
722020	1	● 50	-	50	0.50	0.30	721008	721018	200 350	120
722025	1	● 100	-	100	0.50	0.30	721007	721017	200 350	120 200 300
722030	1	● 200	-	200	0.40	0.20	721007	721017	200 350	200 300
722035	1	● 250	-	250	0.40	0.20	721006	721016	1000	500 1000
722040	1	● 500	-	500	0.30	0.20	721006	721016	1000	500 1000
722045	1	● 1000	-	1000	0.30	0.20	721006	721016	1000	1000
722050	1	● 2000	-	2000	0.30	0.15	721005	721015	5000	
722055	1	● 5000	-	5000	0.30	0.15	721005	721015	5000	



Carousel Stand

Ordering Information

Pipette Stands and Accessories

Cat. No.	Item
725600	Carousel Stand for 6 pipettes
725620	Linear Stand
721259	Pipette Holder for 1 Proline pipette



Linear Stand

Pipette Holder
-Proline Pipettes



■ Pipette Tips

Table of Contents

52 Optifit Tips and SafetySpace™ Filter Tips

54 SafetySpace™ Filter Tips

56 Package Options

62 Tip Compatibility Charts

Optifit Tips and SafetySpace™ Filter Tips

The Perfect Match for Your Pipette



In liquid handling, pipetting results are not dependent on the pipette or the tip alone, but a combination of these and the comfort of the user. Our Optifit Tips and SafetySpace™ Filter Tips are designed and manufactured as a perfect fit for our pipettes, enabling maximal accuracy, precision and ergonomics. The unique Optiload mechanism of our pipettes allows tip attachment and loading with reduced force, but with complete sealing. These tips' universal design also secures compatibility with the pipettes of most other manufacturers.

Manufacturing the tips in our own production facility allows us to maintain the highest quality and purity standards, by selecting the best plastic materials and controlling the manufacturing process from beginning to end. Our quality management system follows not only ISO 9001 and ISO 14001, but also ISO 13485. Tip production also abides by the ISO 14644-1 standard, in order to fulfil ISO Class 8 cleanroom conditions for secured tip purity.



ISO Class 8 Clean Room Manufacturing

Contamination Free Tips

To avoid contamination from human contact, we have automated the entire tip manufacturing process. Pure virgin polypropylene (PP) plastic is automatically fed from silos into moulding machines. Moulding machines and robots located in isolated clean cells load the tips automatically into tip trays and packaging. HEPA filters and higher air pressure are applied for purity within the cell. All Sartorius Single Tray tip racks and Single Refill Packs are individually and automatically packed in air-tight plastic, in order to rule out any danger of contamination.

Additionally, our highly experienced and trained personnel are equipped with specially designed clothing, masks, hair nets and gloves, in order to further diminish contamination risks.

An independent laboratory checks each Single Tray and Refill Pack tip lot for RNase, DNase and endotoxins. Lot-specific purity certificates can be downloaded from www.Sartorius.com select from the navigation bars: -> Liquid Handling -> Tip Quality Certificates.



Fully autoclavable





□ Features and Benefits

Best fit – highest possible accuracy

- **Perfect fitting and sealing with Sartorius pipettes secure the highest possible accuracy and precision**
- Compatible with Optiload feature in Sartorius Picus, eLINE, mLINE and Proline Plus pipettes enabling ergonomic and light tip attachment and ejection
- Colour coding of the tip trays allows easy matching with a suitable colour coded Sartorius pipette
- Compatible with most other pipette makes (see compatibility with other manufacturers' pipettes at www.sartorius.com)

Large Variety of Packaging Options Available:

- Single Tray racks
- Refill Towers
- Single Refill Packs
- Bulk in box


Premium Quality and Purity:

- Strict quality standards, ISO 9001 and ISO 13485, followed from R&D to production and delivery
- Manufactured in ISO 8 classified clean room conditions
- DNase, RNase and endotoxin free manufacturing process:
 - Single Trays and Refill Packs certified pure by lot number
- Pre-sterilised tips are e-beam irradiated
- All tip packages, including individual racks, are lot numbered for full traceability
- The highest quality virgin polypropylene used as raw material



■ SafetySpace™ Filter Tips

No Risk of the Sample Penetrating the Filter



SafetySpace™ Filter Tips leave more space between the sample and the filter than conventional filter tips. This allows pipetting any type of liquid or using any pipetting technique without the risk of the precious sample permeating the filter. In effect, this ensures full sample recovery and thus greater accuracy.

This Feature is Particularly Useful in the Following Applications



- pipetting foaming liquids such as buffers and proteins
- pipetting solvents
- multiple dispensing functions of electronic pipettes
- reverse pipetting

The filter tips are made of virgin polypropylene and feature filter (PE) barriers, in order to prevent aerosol and liquid contamination. This helps protect against the risk of cross-contamination and reduces the pipette's maintenance requirements. The filter is made of polyethylene without "self-sealing" additives to avoid any interference with the sample and results.

SafetySpace™ Filter Tips are available in multiple volumes from 10 µl to 1200 µl and are packed in colour coded Single Tray racks. They are certified DNase, RNase and endotoxin-free and are e-beam pre-sterilised. Lot-specific purity certificates can be downloaded from www.biohit.com.

The SafetySpace™ Filter Tips are Ideal for

- molecular biology
- microbiology
- cell culture applications
- radioactive work



□ Package Options



SafetySpace™ Filter Tips
Pre-sterilised
Single Tray racks



SafetySpace™ Filter Tips

Single Tray Racks

- 10 racks of 96 tips
- Pre-sterilised (e-beam) and DNase, RNase and endotoxin-free certified
- Informative rack labelling improves tip identification and traceability: volume, product number, lot number
- Racks, trays and tips are 100% recyclable.
- Air-tight plastic wrapping around the rack secures purity during transport and storage (the plastic wrapping is regular waste)
- Tip trays are colour coded to indicate the matching, colour coded Sartorius pipette
- Covers a large range of tip volumes from 10 µl to 1200 µl
- Trays and racks are fully autoclavable at 121°C for 20 minutes



Optifit Tips
Single Tray racks



Optifit Tips

Single Tray Racks

- 10 racks of 96 tips
- DNase, RNase and endotoxin-free certified
- Available as e-beam pre-sterilised or non sterilised
- Informative rack labelling improves tip identification and traceability: volume, product number, lot number
- Racks, trays and tips are 100% recyclable.
- Air-tight plastic wrapping around the rack secures purity during transport and storage (the plastic wrapping is regular waste)
- Tip trays are colour coded to indicate the matching, colour coded Sartorius pipette
- Covers a large range of tip volumes from 10 µl to 5 ml
- Fully autoclavable at 121°C for 20 minutes



Refill Towers

Optifit Tips

Refill Towers

- Space-saving with 10x96 tips in one tower.
- Ecological with cardboard packaging, 100% recyclable
- Tip trays to be transferred to Single Tray racks before use
- Trays are colour coded to indicate the matching, colour coded Sartorius pipette
- Covers the most used tip sizes: 10 μ l, 200 μ l and 350 μ l
- Trays and tips are fully autoclavable at 121°C for 20 minutes

Single Refill Packs
Pre-sterilised Optifit Tips

Single Refill Packs

- 10, 15 or 20 trays of 96 tips depending on tip volume
- DNase, RNase, endotoxin-free certified
- Available as e-beam pre-sterilised or non-sterilised
- Individually packed air-tight tip trays for maximum purity with less packaging material compared to racked tips
- Trays to be transferred to Single Tray racks before use
- Container lid with informative tip identification markings including product and lot numbers
- Tip trays are colour coded to indicate the matching, colour coded Sartorius pipette
- Covers a large range of tip volumes from 10 μ l to 1200 μ l
- Trays and tips are fully autoclavable at 121°C for 20 minutes
- Apart from the container lid, material is recyclable








Bulk in a box




Bulk in a Box

- Tips made to the Sartorius quality standard in economical packaging
- Packed in releasable bags in cardboard boxes (1000, 400, 250 or 100 pcs depending on tip volume)
- 100% recyclable
- Covers a large range of tip volumes from 10 μ l to 10 ml
- Tips are fully autoclavable at 121°C for 20 minutes








Ordering Information

Optifit Tips

Volume Range		Length	Packaging	Purity Level RNase, DNase, endotoxin- free	Pre- steri- lised	Qty/Unit	Cat. No.
● 0.1-10 µl		31.5 mm	Single Tray	•		10×96	790010
			Single Tray	•	•	10×96	790011
			Refill Tower			10×96	790012
			Refill Pack	•	•	20×96	790013
			Bulk in Box			1000	790014
● 0.1-10 µl	Extended	46 mm	Single Tray	•		10×96	783210
			Single Tray	•	•	10×96	783211
● 0.5-200 µl		51 mm	Single Tray	•		10×96	790200
			Single Tray	•	•	10×96	790201
			Refill Tower			10×96	790202
			Refill Pack	•	•	15×96	790203
			Bulk in Box			1000	790204
● 5-350 µl		54 mm	Single Tray	•		10×96	790350
			Single Tray	•	•	10×96	790351
			Refill Tower			10×96	790352
			Refill Pack	•	•	15×96	790353
			Bulk in Box			1000	790354
● 10-1000 µl		71.5 mm	Single Tray	•		10×96	791000
			Single Tray	•	•	10×96	791001
			Refill Pack	•		10×96	791002
			Refill Pack	•	•	10×96	791003
			Bulk in Box			1000	791004
			Bulk in Box			400	791005
● 10-1000 µl	Wide bore tip	68.5 mm	Single Tray	•		10×96	791020
			Single Tray	•	•	10×96	791021
			Bulk in Box			1000	791024
● 50-1200 µl		71.5 mm	Single Tray	•		10×96	791200
			Single Tray	•	•	10×96	791201
			Refill Pack	•		10×96	791202
			Refill Pack	•	•	10×96	791203
			Bulk in Box			1000	791204
● 50-1200 µl	Extended	90 mm	Single Tray	•		10×96	791210
			Single Tray	•	•	10×96	791211
			Refill Pack	•		10×96	791212
			Refill Pack	•	•	10×96	791213

Volume Range	Length	Packaging	Purity Level RNase, DNase, endotoxin- free	Pre-steri- lised	Qty/Unit	Cat. No.
● 100-5000 µl	150 mm	Single Tray	•		50	780304
		Single Tray	•	•	50	780305
		Bulk in Box			100	780300
		Bulk in Carton			1000	780308
						
● 1-10 ml	155 mm	Bulk in Box			250	LH-780316
						
0.5-10 ml	Extra long tip for MidiPlus	242 mm	Bulk in Box		100	780310
						

SafetySpace™ Filter Tips

Volume Range	Length	Packaging	Purity Level		Qty/Unit	Cat. No.
			RNase, DNase, endotoxin- free	Pre- sterilised		
● 0.1-10 µl 	32 mm	Single Tray	•	•	10x96	790011F
● 0.5-20 µl 	51 mm	Single Tray	•	•	10x96	790021F
● 2-120 µl 	51 mm	Single Tray	•	•	10x96	790101F
● 5-200 µl 	52.5 mm	Single Tray	•	•	10x96	790201F
● 5-300 µl 	52.5 mm	Single Tray	•	•	10x96	790301F
● 50-1000 µl 	78 mm	Single Tray	•	•	10x96	791001F
● 50-1200 µl 	90 mm	Single Tray	•	•	10x96	791211F

For your guidance the tips are illustrated in the actual size.

Filter tips are not recommended to be used simultaneously with Safe-Cone Filters

Filter Tips with Regular Air Gap

Volume Range	Length	Packaging	Purity Level		Qty/Unit	Cat. No.
			RNase, DNase, endotoxin- free	Pre- sterilised		
● 0.1-10 µl	46 mm	Single Tray	•	•	10x96	783201
						
● 10-500 µl	78 mm	Single Tray	•	•	10x100	783206
						

Empty Tip Boxes for Refill System (Tips and Trays are Not Included)

Item	Tip Type (Non-filtered Tips)	Qty/Unit	Cat. No.
Empty Tip Box for Refill System	10, 200, 350 µl	10	790910
Empty Tip Box for Refill System	1000, 1200 µl	10	790920

Tip Compatibility Charts

Mechanical Pipettes			Optifit Tips, Non-filtered							SafetySpace™ Filter Tips							
			10*	10 Ext	200	350	1000	1000 WB	5000	10 000	10*	10 Ext**	20	120	200	300	500**
	Cat.-No.	µl															
mLINE 1-channel	725010	0.1-3	•	•							•	•					
	725020	0.5-10	•	•							•	•					
	725030	2-20			•								•				
	725050	10-100			•	•							•				
	725060	20-200			•	•								•	•		
	725070	100-1000					•	•									•
	725080	500-5000							•								
	725090	1-10 ml								•							
mLINE 8-channel	725120	0.5-10	•	•							•	•					
	725130	5-100			•	•							•				
	725140	30-300				•									•		
mLINE 12-channel	725220	0.5-10	•	•							•	•					
	725230	5-100			•	•							•				
	725240	30-300				•									•		
Proline Plus 1-channel	728010	0.1-3	•	•							•	•					
	728020	0.5-10	•	•							•	•					
	728030	2-20			•								•				
	728040	5-50			•								•				
	728050	10-100			•	•							•				
	728060	20-200			•	•								•	•		
	728070	100-1000					•	•									•
	728080	500-5000							•								
	728090	1-10 ml								•							
Proline Plus 8-channel	728120	0.5-10	•	•							•	•					
	728130	10-100			•	•							•				
	728140	30-300				•									•		
Proline Plus 12-channel	728220	0.5-10	•	•							•	•					
	728230	10-100			•	•							•				
	728240	30-300				•									•		

* Extended inaccuracy and imprecision with pre-sterilized 10 µl tips

** Filter tips with regular air gap

			Optifit Tips, Non-filtered							SafetySpace™ Filter Tips								
			10*	10 Ext	200	350	1000	1000 WB	5000	10 000	10*	10 Ext**	20	120	200	300	500**	1000
	Cat.-No.	µl																
Mechanical Pipettes Proline Plus Fixed Volume 1-channel	728515	5	•	•														
	728520	10	•	•														
	728530	20			•							•						
	728535	25			•								•					
	728545	50			•								•					
	728550	100			•	•							•					
	728560	200			•	•								•	•			
	728565	250					•	•									•	•
	728567	500					•	•									•	•
	728570	1000					•	•										•
	728575	2000																•
	728580	5000																•
	728590	10 ml																•
Proline 1-channel	720005	0.1-2.5	•	•														
	720000	0.5-10	•	•	•													
	720080	2-20			•							•	•					
	720020	5-50			•	•							•					
	720050	10-100			•	•							•	•	•			
	720070	20-200			•	•								•	•			
	720060	100-1000					•	•										•
	720110	1000-5000																•
Proline 8-channel	720210	0.5-10	•	•														
	720220	5-50			•	•							•	•	•			
	720240	50-300				•									•			
Proline 12-channel	720310	0.5-10	•	•														
	720320	5-50			•	•							•	•	•			
	720340	50-300				•									•			
Proline Fixed Volume 1-channel	722001	5	•	•	•								•	•	•			
	722004	10	•	•	•								•	•	•			
	722010	20			•	•							•	•				
	722015	25			•	•							•					
	722020	50			•	•							•					
	722025	100			•	•							•	•	•			
	722030	200			•	•								•	•			
	722035	250					•	•									•	•
	722040	500					•	•									•	•
	722045	1000					•	•										•
	722050	2000																•
	722055	5000																•

* Extended inaccuracy and imprecision with pre-sterilized 10 µl tips

** Filter tips with regular air gap

Electronic Pipettes			Optifit Tips, Non-filtered										SafetySpace™ Filter Tips							
			10*	10 Ext	200	350	1000	1000 WB	1200	1200 Ext	5000	10 000	10*	10 Ext**	20	120	200	300	500**	1000
	Cat. No.	µl																		
Picus 1-channel	735021	0.2-10	•	•											•	•				
	735041	5-120			•	•											•	•		
	735061	10-300				•												•		
	735081	50-1000					•	•												•
	735101	100-5000									•									
	735111	500-10000												•						
Picus 8-channel	735321	0.2-10	•	•											•	•				
	735341	5-120			•	•										•				
	735361	10-300				•											•			
	735391	50-1200							•	•										•
Picus 12-channel	735421	0.2-10	•	•											•	•				
	735441	5-120			•	•										•				
	735461	10-300				•											•			
	735491	50-1200							•	•										•
eLINE 1-channel	73001X	0.1-5	•	•											•	•				
	73002X	0.2-10	•	•											•	•				
	73004X	5-120			•	•										•				
	73006X	10-300				•											•			
	73008X	50-1000					•	•												•
	73010X	100-5000									•									
eLINE 8-channel	73032X	0.2-10	•	•											•	•				
	73034X	5-120			•	•										•				
	73036X	10-300				•											•			
	73039X	50-1200							•	•										•
eLINE 12-channel	73042X	0.2-10	•	•											•	•				
	73044X	5-120			•	•										•				
	73046X	10-300				•											•			
	73049X	50-1200							•	•										•

* Extended inaccuracy and imprecision with pre-sterilized 10 µl tips

** Filter tips with regular air gap





■ Dispensers and Dispenser Tips

Table of Contents

68 eLINE® Lite and Pro Dispenser

70 Dispenser Tips

72 Mechanical Stepper

■ eLINE® Lite and Pro Dispensers

A Single Innovation for Repetitive Dispensing

eLINE Pro and Light electronic dispensers are specifically designed for convenient and ergonomic repetitive liquid dispensing. Together with dispenser tips, these dispensers operate according to the positive displacement principle. They are therefore excellent tools for dispensing viscous or infectious liquids, or liquids with high vapour pressure. Their unique TipGuide feature selects and displays the correct tip size saving time, and eliminating the risk of choosing the incorrect tip. eLINE Dispensers also include the innovative electronic tip ejector feature, familiar from the Sartorius range of electronic pipettes, for reducing strain on the user's hand.

- Eases the work load in long pipetting series
- High accuracy in dispensing viscous or infectious liquids or liquids with high vapour pressure
- Contamination-free dispensing based on positive displacement principle
- Quick and error-free tip selection with unique, built-in TipGuide, which suggests the correct tip size automatically
- Large volume range from 1 µl to 50 ml
- Compatible with Sartorius electronic pipette charging stands
- Compatible with Sartorius Dispenser Tips

TipGuide selects and displays the suitable tip size

- No time-consuming volume calculations
- Just set the desired volume and number of dispensings

Direct charging option while dispensing

Clear LCD-display with charging indicator

Easy-to-use keypad enables fast programming

Unique electronic one-touch tip ejection

Electronic start button for all pipetting operations

Attractive ergonomic design

Both left and right-handed operation

Easy and light tip insertion

Dispenser Tips

- Allow contamination free dispensing





Operating modes and functions

	eLINE Pro	eLINE Lite
Multiple Dispensing with TipGuide	✓	✓
Multiple Dispensing	✓	
Pipetting	✓	
Diluting	✓	
Sequential Dispensing	✓	
Automatic Multiple Dispensing	✓	
Multi-Aspirating	✓	
Custom mode (CST)		
– GL (select service intervals)	✓	
– Sr (select the desired tip range)		
– SET (returns default setting)		
5 speeds for aspirating and dispensing	✓	✓



Ordering Information

eLINE® Dispenser

Cat. No.	Item
73070x	eLINE Pro Dispenser
73080x	eLINE Lite Dispenser
730981	Charging Stand
730991	Charging Carousel for 4 pipettes
731001	eLINE replacement battery

x:=0 without AC-adaptor, x=1 with AC-adaptor (EU, UK, US|JPN, AUS and CHN plugs)

Dispenser Tips



Dispenser tips:

- Function according to positive displacement principle
- Contamination-free dispensing – no aerosols are formed
- Made of virgin polypropylene (tip) and polyethylene (plunger)
- 9 different tip sizes from 0.1 to 50 ml
- DNA, RNase, ATP and endotoxin-free
- Pre-sterilised tips are individually wrapped
- Suitable for use together with eLINE Pro and Lite Dispensers and Sartorius Mechanical Stepper, as well as most other mechanical steppers (Ripette (Ritter), HandyStep (Brand), Minilab 100/101 (HTL), EasyStep (Tomos))

Ordering Information

Dispenser Tip

Cat. No.	Volume (ml)	Purity Level	Qty/Unit
792017	0.1	Non sterile	100
792026		Pre-sterilised	100
792018	0.2	Non sterile	100
792027		Pre-sterilised	100
792019	0.5	Non sterile	100
792028		Pre-sterilised	100
792020	1.0	Non sterile	100
792029		Pre-sterilised	100
792021	2.5	Non sterile	100
792030		Pre-sterilised	100
792022	5.0	Non sterile	100
792031		Pre-sterilised	100
792023	10	Non sterile	100
792032		Pre-sterilised	100
792024	25	Non sterile*	25
792033		Pre-sterilised	25
792025	50	Non sterile*	25
792034		Pre-sterilised	25
792038 Starter Kit	20 x 0.5 ml 20 x 1.0 ml 20 x 2.5 ml 20 x 5.0 ml 20 x 10.0 ml	Non sterile*	100
792036	Adapter for 25 and 50 ml tips, 3 pcs./unit		

* Includes 1 adapter/unit



□ Performance Specifications

eLINE® Pro and Lite Dispenser and Dispenser Tip System

Tip Volume (ml)	Volume Min Max (µl)	Increment (µl)	Step Size Max (µl)	Test Volume (µl)	Inacc. (%)	Impr. (%)	Number of Dispensings Min Max
0.1	1 100	0.2	100	100	1.00	0.50	1 100
				10	1.00	2.00	
0.2	2 200	0.4	200	200	1.00	0.50	1 100
				20	1.70	2.50	
0.5	5 500	1	500	500	0.80	0.60	1 100
				50	1.00	3.00	
1.0	10 1000	2	1000	1000	0.50	0.20	1 100
				100	0.50	0.60	
2.5	25 2500	5	2500	2500	0.60	0.20	1 100
				250	0.60	0.40	
5.0	50 5000	10	5000	5000	0.60	0.15	1 100
				500	0.90	0.50	
10.0	100 10000	20	10000	10000	0.40	0.20	1 100
				1000	0.40	1.00	
25.0	500 25000	50	25000	25000	0.50	0.15	1 50
				2500	0.50	0.50	
50.0	1000 50000	100	50000	50000	0.30	0.15	1 50
				5000	0.30	0.70	

The performance specifications apply to the eLINE Dispenser and Sartorius Dispenser Tip System. Due to Sartorius' continuous product development, specifications may change without prior notice.

■ Mechanical Stepper



Sartorius Mechanical Stepper is an easy-to-use positive displacement dispenser, which allows rapid multi-dispensing of pre-set volumes up to 48 times in succession, without refilling.

The unit comes complete with an adapter for 25 ml and 50 ml tips.

- Lightweight (105 g) yet robust construction
- Maintenance-free design
- Dispensing volumes from 1 μ l to 50 000 μ l (min. volume with Sartorius Dispenser Tip is 2 μ l)
- Ideal for dispensing aqueous and viscous liquids
- Contamination-free dispensing based on positive displacement principle
- Its ergonomic design and positioning of the dosage button in the upper part of the device allows single-handed volume selection, loading and dispensing
- 48 dispensing steps at intervals of 1 second without the need to refill thus saving 90% of the working time required for the usual pipetting technique
- Tested according to ISO 8655-1 standard

□ Performance Specifications

Mechanical Stepper and Dispenser Tips

Tip Volume (ml)	Dosage Volume (μ l)					Impr. (%)	Inacc. (%)
	1	2	3	4	5		
Adjustment on dispenser	1	2	3	4	5		
Number of steps	48	23	15	11	8		
0.10	2	4	6	8	10	< 1.6	< \pm 3.0
0.20	4	8	12	16	20	< 1.3	< \pm 2.0
0.50	10	20	30	40	50	< 0.5	< \pm 0.8
1.00	20	40	60	80	100	< 0.9	< \pm 0.9
2.50	50	100	150	200	250	< 0.4	< \pm 0.8
5.00	100	200	300	400	500	< 0.3	< \pm 0.4
10.00	200	400	600	800	1000	< 0.5	< \pm 0.6
25.00	500	1000	1500	2000	2500	< 0.3	< \pm 0.2
50.00	1000	2000	3000	4000	5000	< 0.2	< \pm 0.2

□ Ordering Information

Mechanical Stepper

Cat. No.	Item
725700	Mechanical Stepper with adapter for 25 ml and 50 ml dispenser tips



■ Stands

Table of Contents

76 Pipette Stands

■ Pipette Stands



Charging Carousel Stand for 4

When the pipette is not in use, it should be stored in an upright position in order to avoid contamination from work surfaces. Sartorius provides stands for all of its pipettes. It is recommended that electronic pipettes be stored and charged on a charging stand whenever they are not in use. In this way, their batteries always remain charged for when work begins.

The Linear Stand is designed for all Sartorius mechanical and electronic pipettes, particularly for mLINE, Proline Plus and Proline mechanical devices. This stand is also compatible with a wide range of other pipette makes.

The simplest of all are the pipette holders which are attached to the front edge of a shelf. These are suitable for mechanical pipettes.



Charging Stand for one

Compact carousel stands are ideal for saving bench space. There is one for mechanical pipettes, and a charging carousel stand for electronic pipettes.



Linear Stand (non-charging)





Non-charging Carousel Stand

□ Ordering Information

Pipette Stands

Cat. No.	Item
730981	Charging Stand for one electronic pipette*
730991	Charging Carousel for 4 electronic pipettes*
725620	Linear Stand for all Sartorius pipettes
725600	Carousel Stand for 6 mechanical pipettes
725610	Holder for one mLINE Proline Plus pipette
721259	Holder for Proline pipette

* Includes Universal AC-adaptor (with Euro, US | Jpn, UK and China plugs)

Pipette Holder
-mLine/Proline Plus PipettesPipette Holder
-Proline Pipettes



■ Accessories

Table of Contents

80 Safe-Cone Filters

82 Elbow Pad

83 Reagent Vessel

83 Cooling Rack

■ Safe-Cone Filters



Why Should You Use Safe-Cone Filters?

These unique and replaceable polyethylene (PE) filters act as a final barrier to prevent any fluids and liquid vapours from reaching the internal components of the pipette.

- Protect the pipette and sample from contamination
- Prolong the pipette's lifetime
- Reduce maintenance intervals
- Cost-effective

When Should You Use them?

The ultimate pipette protectors are available in two types:

- **Plus Filter**
For more demanding applications such as radioactive work, cell culture, bacterial and virological work and molecular biology.
- **Standard Filter**
For general applications. Can be used in same type of work as the Plus filter is recommended for, but needs to be changed more frequently.

How Often Should You Change?

The interval of changing the filter depends completely on the application and the sample. However, according to studies the filter is recommended to be changed daily (after 50 to 250 pipettings) and immediately in case of over-aspiration.


How to Change?

To ensure the safety of the user, forceps should be used to avoid touching the contaminated filters by hand. The mLINE also features a built-in filter ejector. In addition, clean the tip cone with ethanol (70%) prior to the assembly of a new filter.

Contamination in Pipette Barrel

Number of Pipettings	No Filter	Safe-Cone Filter
50	-	-
100	+	-
250	++	-
500	+++	-

+: DNA (50 µl plasmid DNA 120 µg/ml) contamination in pipette barrel.



Pipette contamination in pipette barrel when pipetting liquid culture of bacteria *Micrococcus Luteus*.

Ordering Information

Safe-Cone Filters

Cat. No.	Item	Qty/Unit
721008	Standard Ø 2.51 mm PE	50
721007	Standard Ø 3.15 mm PE	50
721006	Standard Ø 5.33 mm PE	50
721005	Standard Ø 6.73 mm PE	50
721014	Standard Ø 1.83 mm PE	50
721018	Plus Ø 2.51 mm PE	50
721017	Plus Ø 3.15 mm PE	50
721016	Plus Ø 5.33 mm PE	50
721015	Plus Ø 6.73 mm PE	50

PE = polyethylene

■ Elbow Pad



The Elbow Pad helps you to feel more comfortable while pipetting. The visco-elastic material of the pad relieves contact stress, pain and discomfort under your elbow.



The Pad is Ideal for

- long pipetting series
- work requiring high concentration, e.g. micro plate work
- any work where you need a cushion beneath your elbow or wrist



□ Features and Benefits

- Improves pipetting ergonomics
- Forms according to any elbow size or shape
- Coating is pleasant to the skin
- The compact size takes up little desk space
- Very durable
- Easy to clean with washing up liquid, or ethanol (70%)
- Non-autoclavable

□ Ordering Information

Elbow pad

Cat. No.	Item	Qty
723103	Elbow Pad	1

■ Reagent Vessel



Made from polypropylene, the autoclavable and durable reagent vessel is chemically resistant to all common reagents.

□ Ordering Information

Reagent Vessel

Cat. No.	Item	Qty
783500	Reagent Vessel (capacity 120 ml)	16

■ Cooling Rack



The Cooling Rack keeps the reagents cool during pipetting, and is especially suitable for applications in molecular biology. It is compatible with 1.5 ml and 2.0 ml conical and cylindrical microcentrifuge tubes.

□ Ordering Information

Cooling Rack

Cat. No.	Item	Qty
723102	Cooling Rack	1



■ Maxi-volume Liquid Handling

Table of Contents

86	Midi Plus™
88	Biofiller
90	Proline® Prospenser
91	Prospenser
92	Biotrate Digital Burette

■ Midi Plus™

Excellent Performance and Ergonomics

The Midi Plus is a lightweight electronic cordless pipetting controller, which allows aspiration from bottles and tubes, without the arm and hand elevations required in the case of serological or volumetric pipettes. It fits all commonly used glass or plastic pipettes, but can also be used with Sartorius 5 ml and 10 ml disposable tips. The speed can be fine-tuned by applying varying finger pressure to the operating buttons.

The Midi Plus is ideal, for example, in microbiological work: dispensing into a culture media dish can be performed carefully, drop by drop, without breaking the fine surface of the medium.

- Precise finger tip control of variable speed ranges
- Hydrophobic autoclavable filter prevents over-aspiration
- Convenient fold-out bench stand supports the unit and pipette when not in use*
- Linear speed control
- Low battery warning
- CE marked
- 2-year warranty

Hydrophobic autoclavable filter prevents over-aspiration

Autoclavable pipette adapter set

Fits all commonly used 1-100 ml glass or plastic pipettes

Operating buttons



Stepless speed control

- Allows the aspiration and dispensing speed to be adjusted easily and precisely to suit either large or small volume pipettes.

Fold-out stand*

- Prevents pipette contamination
- Saves space
- Allows the unit to be rested on a table with a pipette attached



*Not available in units sold in the United States

□ Features

Pipette types	Plastic or glass 1– 100 ml Pasteur pipettes 5 ml and 10 ml Sartorius pipette tips
Rechargeable during use	Yes
Speed control	Stepless adjustable control
Gravity dispensing	Yes
Stand	Attached support
Weight	207 g
Low battery indicator	Yes
Autoclavable parts	Nose cone, pipette holder and filter

□ Ordering Information

Cat. No.	Item	Qty/Unit
710931	Midi Plus Pipetting Controller with Universal AC-adaptor ¹	1
712912	Replacement filter, 0.45 µm, autoclavable	5
712913	Replacement filter, sterile, 0.45 µm	1
711015 ²	Adapter set for 5 ml tip, autoclavable	1
711016 ²	Adapter set for 10 ml tip, autoclavable	1
780300	Optifit Tip 5 ml	100
780310	Midi Plus Tip 10 ml	100

1) Universal AC-adaptor (incl. Euro, US | Jpn, UK and China plugs)

2) Adapter set including the nose cone and the silicon adapter



■ Biofiller

For Smooth Manual Pipetting

Biofiller is a uniquely designed, lightweight and easy-to-use pipette filler with an ultra-squeezable bulb to provide smooth, manual pipetting control when both aspirating and dispensing, using 1 ml to 100 ml pipettes.

Simply squeeze the large silicone bulb, and the thumb lever controls both the aspirating and dispensing modes, with a button for blowing out residual contents if required.

- Comfortable and simple to use
- Precise pipetting control
- Robust and lightweight
- Compatible with blow-out pipettes
- Uses integral 0.45 µm, replaceable membrane filter to ensure liquid is not accidentally drawn into the unit
- Autoclavable silicone pipette holder





□ Ordering Information

Cat. No.	Item	Qty/Unit
723039	Biofiller pipette filler with integral 0.45 µm filter	1
721963	Replacement silicone pipette holder	1
721962	Replacement filter, 0.45 µm	5
721965	Replacement filter, 0.20 µm	5
721966	Replacement filter, sterile, 0.45 µm	1

Proline® Prospenser

Easy-to-use Bottle-top Dispenser

Proline Prospenser has been designed for trouble-free and reliable bottle-top dispensing of liquids, including concentrated acids, bases, saline solutions, as well as many organic solvents.

- Dispensing directly from the supply bottle
- Fully autoclavable at 121 °C, 2 bar, 20 minutes
- Accuracy of delivery within $\pm 0.5\%$
- Chemically resistant fluid pathway
- Anti-drip closure cap included
- Easy-to-use volume adjustment for reproducible dispensings
- Easy to dismantle for cleaning and maintenance
- Wide range of adapters included to fit the most common bottle sizes
- Optional flexible dispense tube extension (max length 800 mm) with safety handle enables fast and precise dispensing even into narrow tubes
- Each unit is supplied with performance certificate and tested according to ISO 8655

Volume setting easy to adjust

Borosilicate glass barrel protected with a transparent polypropylene sleeve

Anti-drip closure cap

Wide range of bottle adapters included

Chemically resistant liquid pathway



Optional flexible dispense tube extension
- Coiled
- Length 800 mm



Fully autoclavable

Prospenser

Bottle-top Dispenser with Anti-drip Valve

The Prospenser bottle-top dispenser delivers safe and precise liquid dispensing, including of strong acids, bases and solvents. Special features include an easy-to-adjust calibration mechanism and precision valve for enhanced accuracy and usability.

- Dispensing directly from the supply bottle
- Easy-to-adjust calibration mechanism
- Fully autoclavable at 121°C, 2 bar, 20 minutes
- Accuracy of delivery within $\pm 0.3\%$
- Chemically resistant fluid pathway
- Anti-drip precision valve mechanism ensures easy priming and minimum waste with no leakage back into the reservoir
- Easy-to-use volume adjustment for reproducible dispensing
- Unlike other bottle-top dispensers, Prospenser's glass barrel can be disassembled from the pedestal for thorough cleaning
- Wide range of adaptors included to fit the most common bottle sizes
- Optional dispense tube extension allows fast and safe dispensing even into narrow tubes
- Manufactured to ISO9002 standards, each unit is supplied with an individual calibration certificate



Easy-to-adjust calibration mechanism.



Borosilicate glass barrel protected with a transparent polypropylene sleeve

Easily removable PTFE piston for cleaning and smooth action

Chemically resistant liquid pathway

Precision valve mechanism
- Ensures that the Prospenser stays fully primed all day

Bubble free dispensing

Anti-drip tap valve

Set of adaptors available



Optional dispense tube extension
- Allows fast and safe dispensing into narrow tubes



Fully autoclavable

■ Biotrate Digital Burettes

Ranges of 0–30 ml and 0–50 ml

The streamlined Biotrate digital burette and dispenser delivers accurate, precise and convenient bottle-top titration, as well as optimum operator safety. Due to its life-time battery concept, there is no need for an electrical outlet. This feature makes Biotrate ideal for both laboratory and field-based analyses, where it can be easily moved from one place to another.

- Chemically resistant and autoclavable liquid-path construction is made of borosilicate glass, PTFE, PVDF, FEP and ceramic components
- Based on positive displacement principle for maximum accuracy
- An automatic low battery indication is shown on the large easy-to-read LCD display
- Simple for user to calibrate and restore factory calibration
- With a zero reset feature, it is easy to move from one titration to another



Rotating head and clear numbers on the display increase working reliability

Chemically resistant and autoclavable liquid pathway

Simple user calibration and factory reset for ISO 9000 and GLP requirements

Visible borosilicate glass barrel protected with a transparent polypropylene sleeve allows inspection of liquid before dispensing

Bubble-free dispensing

Anti-drip tap



□ Ordering Information and Performance Specifications

Proline® Prospenser

Cat. No.	Item	Increment	Max Volume	Inacc. (%)	Impr. (%)
723045	Proline Prospenser 0.5 – 5 ml (with 25, 28, 32, 38 and 40 mm adaptors)	0.1 ml	5 ml	0.5	0.1
723046	Proline Prospenser 1 – 10 ml (with 25, 28, 32, 38 and 40 mm adaptors)	0.2 ml	10 ml	0.5	0.1
723047	Proline Prospenser 2.5 – 25 ml (with 32, 38 and 40 mm adaptors)	0.5 ml	25 ml	0.5	0.1
723048	Proline Prospenser 5 – 50 ml (with 32, 38 and 40 mm adaptors)	1.0 ml	50 ml	0.5	0.1
721633	Flexible tube extension for 5 and 10 ml Proline Prospensers				
721634	Flexible tube extension for 25 and 50 ml Proline Prospensers				

Prospenser

Cat. No.	Item	Increment	Max Volume	Inacc. (%)	Impr. (%)
723049	Prospenser 0.01 – 2.5 ml (with 38, 40 and 45 mm adaptors)	0.05 ml	2.5 ml	0.3	0.1
723050	Prospenser 0.1 – 5 ml (with 38, 40 and 45 mm adaptors)	0.1 ml	5 ml	0.3	0.1
723051	Prospenser 0.2 – 10 ml (with 38, 40 and 45 mm adaptors)	0.2 ml	10 ml	0.3	0.1
723052	Prospenser 1 – 30 ml (with 38, 40 and 45 mm adaptors)	1.0 ml	30 ml	0.3	0.1
723053	Prospenser 1 – 50 ml (with 38, 40 and 45mm adaptors)	1.0 ml	50 ml	0.3	0.1
721998	Dispense tube extension				

Biotrate Digital Burettes 0 – 30 ml and 0 – 50 ml

Cat. No.	Item	Increment	Inacc. (%)	Impr. (%)
723054	Biotrate 0 – 30 ml (with 33, 38 and 45 mm adaptors)	0.01 ml	0.2	0.1
723055	Biotrate 0 – 50 ml (with 33, 38 and 45 mm adaptors)	0.01 ml	0.2	0.1
721998	Expandable delivery jet			



■ Pipetting Academy™


Table of Contents

96 Pipetting Academy™

98 Pipetting Recommendations

■ Pipetting Academy™

Training for Better Performance, Ergonomics and Safety



Seminars include both theory and practice!

Videos and animations support learning.

Are you concerned about the results of your work due to poor pipetting practices or RSI (Repetitive Strain Injury)?

Have you considered that the pain in your hand or arm may be related to the instruments or techniques you use?


Do you know which pipetting technique to use with different types of liquids?

Do your results vary between users?

The Pipetting Academy™ seminar offers a comprehensive tutorial package, developed to answer these questions together with you. During the seminars, you will learn to recognise pipetting-related risk factors and increase your knowledge of ergonomics, safety and pipetting techniques, in order to avoid these risks in your daily work.

What Will You Learn?

- Get hands-on training in pipetting techniques that will help you and your co-workers to obtain more accurate and precise results in the lab
- Gain a better understanding of the influence of pipetting techniques and environmental factors on testing results
- Be guided through the essentials of laboratory ergonomics
- Gain a better understanding of the ergonomic risks in the laboratory environment and liquid handling in particular
- Learn how to avoid these risks by choosing the most appropriate working postures, liquid handling devices and accessories
- Appreciate how you can help make savings in both direct and indirect costs due to bad ergonomics
- Become able to instruct your employees on all of these issues, making work more efficient and enjoyable.



Build Your Own Seminar

Pipetting Academy offers you various seminars for different purposes.

You May Choose from

Ergonomics

Learn about the optimal posture for pipetting and become familiar with tools that can help you work ergonomically and efficiently. Understand the risks and learn about the solutions.

Pipetting Techniques

Master your working tool. Handle the pipette correctly. Be guided through the many techniques of which your pipette is capable.



How to Sign Up for Seminars?

- To sign up for the seminar, contact your local Sartorius representative
- The seminar will be held in the location most suitable to you and your colleagues
- The trainer will be certified to hold Pipetting Academy seminars
- Each participant will receive a certificate of participation after the seminar

Gain access to educational material, videos and animations

Once you have signed up and participated in the seminar, you will automatically gain access to educational videos, animations and presentations on ergonomics, pipetting techniques and calibration.



Pipette Calibration

Learn how to calibrate your pipette correctly to obtain the most accurate pipetting results.


ISO 8655 Quality Standards

All you need to know about the standards you face in your everyday work.


- Through these animations, you will be guided step by step in the correct handling of the pipette and through reverse pipetting, diluting and all other pipetting modes, to make your work easier and more efficient
- A Certified Professional Ergonomist will guide you through the essentials of ergonomics
- Presentations with explanations and illustrations on calibration and quality standards, are also available

See You at the Pipetting Academy™!

■ Pipetting Recommendations



Hold the pipette in a vertical position when pipetting



Pre-rinsing the tip is recommended if the temperature of the liquid is the same as the ambient temperature



Avoid contamination with Safe-Cone Filters



mLINE volume lock prevents volume changes during pipetting

Preparations Before Pipetting

- Use the tip specified by the manufacturer.
- Ensure that the pipette and the tip have been tested according to ISO 8655 and the tip is seated correctly.
- Make sure pipettes have been correctly calibrated.
- Check that the pipette, tip and liquid are all at the same temperature.
- When pipetting liquids with temperatures different to the ambient temperature, do not pre-rinse the tip. Change the tip after each pipetting.
- Ensure that any fluid viscosity variations have been accounted for and the correct technique is employed, i.e. reverse pipetting.
- If handling infectious or radioactive agents make sure appropriate shielding and other precautions protect the operator.
- Use Safe-Cone Filter in the tip cone whenever possible.

While Pipetting

- Hold the pipette in a vertical position. Tilting the pipette at an angle causes a volume greater than the set volume of liquid to enter the tip.
- In most cases, pre-rinsing of the tip is recommended, to achieve accurate results. Do not pre-rinse the tip, if the temperature of the liquid is different to the ambient temperature.
- When aspirating fluid, the pipette tip should normally be immersed to a depth of 2-3 mm.
- When using a mechanical pipette, operate the piston with a smooth and consistent thumb action, for repeated results without foaming or bubbles.
- You should pipette against the inside wall of the receiving vessel. Remove the tip by drawing it upwards against the inside wall.
- Ensure that the pipette blow out action is fully activated.
- Ensure that the volume is still set at the required position. A pipette with a volume locking mechanism is recommended, in order to avoid accidental volume change during pipetting.
- Avoid leaving the pipette on its side with liquid in the tip, which may seep back into the mechanism.



Charging while pipetting is possible with Sartorius electronic pipettes



Load the tip onto the pipette carefully and take advantage of the Optiload tip loading mechanism



Clean the pipette before sending it to service

Other Precautions

- Store the pipette on a stand when not in use – see our product pages for more information on pipette stands. Electronic pipettes should be returned to their charging stands.
- Avoid dropping the pipette or allowing contact with dirt or grease.
- Change the Safe-Cone Filter regularly (recommendation after 50 to 250 pipetting cycles), and in every case of over-aspiration.
- Never strike the tip cone against the tip tray when loading the tip, as this can damage the pipette.
- Avoid exposing the unit to extreme temperature changes, humidity and dust (operating temperature from 15°C to 40°C).
- Service the pipette regularly.
- Clean the pipette thoroughly before sending it in for service. Decontaminate the pipette with 70% ethanol. Notify the service personnel of the purpose for which the instrument has been used. Postal services may refuse to deliver instruments used for hazardous materials. Make sure that a qualified person services the pipette.



■ Pipette Calibration and Maintenance Services

Table of Contents

- 102 Pipette Calibration and Maintenance Services
- 108 Quanta – Pipette Service and Calibration Software
- 112 Pipette Decontamination Procedure
- 113 Autoclaving Instructions
- 114 Troubleshooting Guide

■ Pipette Calibration and Maintenance Services



Why is Calibration and Maintenance needed?


All pipette manufacturers recommend that regular maintenance and calibration is applied to maintain reliable pipetting results and to maximize the lifetime of the pipette.

Pipette calibration is a fundamental part of Good Laboratory Practice (GLP) and quality systems and must be considered a vital part of any laboratory regime where precise volumes of fluid need transferring or diluting. Pipette performance is measured as accuracy and precision or how close the dispensed volumes are to the target and how close the results are to each other.

Sartorius Service Centers

It is vitally important to Sartorius that our customers receive world class service and support, from the first phone call to the moment the engineer leaves, with the customer's equipment in perfect working condition.

Sartorius provides a global network of service centers for calibration of all makes and models of pipettes, burettes, bottle top dispensers and other liquid handling instruments. With over 20 years of experience in pipettes and liquid handling instrument services, the global organization provides world class services tailored to customers' needs, simultaneously fulfilling standards and regulatory body requirements. With the Sartorius concept of "all makes and models liquid handling services" you can be certain that your instruments are handled according to the international ISO 8655 standard defining the use and calibration of pipettes and other liquid handling instruments. Sartorius provides reliability of operation, reduced instrument downtime and confidence, so that your work is according to the strictest requirements.





□ Maintenance, Repair and Calibration Services

Services Provided

- Maintenance and calibration for all makes and models of pipettes and other liquid handling instruments, including cleaning, relubrication and adjustment of instruments
- Pipette calibration done according to ISO 8655 standard by our network of ISO 17025 accredited pipette calibration laboratories around the world
- Pipette maintenance and calibration services available either as mail-in service at service centers or as on-site service at the customer's laboratory
- Repair services and spare parts for all makes and models of pipettes available
- Extended Warranty for an additional 12 months

Benefits

- Maximum reliability and lifetime with regular maintenance and original spare parts
- Confidence in the quality of your work provided by regular calibration
- All documentation needs met by our detailed service and calibration documentation, including ISO 17025 accredited calibration certificates with measurement results traceable to international standards
- Efficient and time saving services provided either at your site or as a mail-in service, which includes courier collection and delivery
- Affordable, planned warranty cover of your equipment



□ Maintenance

Maintenance consists of the cleaning and lubrication of the piston and tip cone(s) of the pipette, and the replacement of parts prone to wear, in order to guarantee a long lifetime and consistently accurate and precise results.

As with all mechanical devices, pipettes may need repairs, and it is usually worth the expense to maintain and repair your pipettes rather than to dispose of them. This is also an ecologically sustainable choice. Sartorius provides you with spare parts and repairs for any make and model. If the pipettes are beyond economical repair, then we offer replacements.

Calibration and maintenance can be done either at a dedicated calibration laboratory, or at the customer site. Laboratory conditions can provide a smaller uncertainty of measurement in calibration, but the quickest turnaround time is often achieved by performing maintenance and calibration on-site.





□ Accreditation and Calibration Standards

Accreditation and Quality Systems

Good Laboratory Practices (GLP), accredited processes and quality systems require traceability of measurement and compliance with a multitude of standards. Sartorius operates a global network of accredited calibration laboratories, providing services around the world.

What is ISO 17025?

The ISO 17025 standard specifies the general requirements for the competence of calibration laboratories to carry out tests and calibrations. A pipette calibration laboratory with an ISO 17025 accreditation provides the highest level of reliability and confidence in pipette calibration, proven with a Measurement Uncertainty Certificate provided with each calibration, when required. Gaining and maintaining an ISO 17025 accreditation requires extensive reviewing and development of personnel, procedures and facilities, audited annually by the accreditation body.

Sartorius provides ISO 17025 accredited pipette calibration services in various countries, for example DAkkS in Germany, Cofrac in France, UKAS in the United Kingdom, A-Class in the USA, JCSS in Japan, and FINAS in Finland. All European

accreditation bodies are members of the European Co-operation for Accreditation (EA) and signatories of the Multilateral Recognition Agreements (MLA) ensuring that accredited pipette calibration certificates are accepted around the world.

Ask your local Sartorius representative for more information on accredited pipette calibration services available for you.

What is ISO 8655?

The ISO 8655 standard specifies the requirements for piston operated volumetric apparatus (pipettes) and pipette calibration laboratories, providing detailed requirements for procedures and equipment used in pipette calibration. Pipettes calibration, according to ISO 8655, is done in a carefully controlled environment, with no drafts or vibrations, using repeatable and reliable measurement technology. Sartorius balances meet such specifications and calibration is always done according to the strictest, repeatable procedures.

Making sure your pipettes are calibrated according to ISO 8655, Sartorius provides the best confidence in the reliability of results. Sartorius pipette calibration is always done according to ISO 8655.





In Conclusion:

Adopting a regular calibration and maintenance routine for your pipette has the following benefits:

Confidence

Your pipettes are operating correctly with the accuracy and precision you know.

Reliability

With maintenance and calibration and operational qualification you are able to trust the instrument's operational capability.

Efficiency and effectiveness

With properly working pipettes you can work uninterrupted and be more efficient.

Sustainability

Extending the lifespan of your pipette reduces waste, and consequently is an environmentally considerate choice.

Frequently Asked Questions

Question: What makes and models of pipettes does Sartorius service?

Answer: Sartorius provides calibration and maintenance services for all manufacturers' pipettes, the most common being Sartorius (*Biohit family*), Gilson, Eppendorf, and Rainin. With over 20 years of experience in the field we have the skill and ability to service any pipette.

Question: What is pipette calibration?

Answer: It is when a pipette is tested for inaccuracy, imprecision and linearity to determine the uncertainties of measurement for each volume calibrated. Test volumes and data points are selected according to the customer requirements.

Question: Can I have my pipettes calibrated in my laboratory, as I need them every day?

Answer: Yes. Sartorius offers calibration and maintenance services both at our workshops and at your laboratory site. Ask for on-site services, especially if you need your pipettes every day.

Question: Why should I get a Service Contract?

Answer: Often the most economical and easiest way to handle your calibration and maintenance needs, is to enter into a contract with a service provider. This saves time, and is usually the most economical choice. Sartorius service technicians are experts in pipettes and also offer training in the use of instruments. They can also offer replacements for worn out pipettes.

Question: How often should my pipettes be calibrated?

Answer: The ISO 8655 standard states that pipettes should be calibrated regularly, at least annually, or more often, for example every 3-6 months. It depends on the accuracy requirements of your work as well as the nature of liquids used - corrosive or volatile liquids create more wear and tear on your pipettes.



■ Quanta – Pipette Service and Calibration Software




Instruments and equipment used within a laboratory must be routinely tested, validated and calibrated by laboratory personnel or an external technically competent and authorised calibration company.

Quanta Range of Software

Sartorius provides anyone who needs to check and calibrate pipettes, from end users to pipette service companies, with special software to ensure the ongoing performance of their pipettes. Quanta schedules pipette calibrations, controls and captures data from balances automatically, performs all required calculations and produces reports such as calibration reports, service reports, overdue pipettes reports etc. Quanta software is suitable for both single and multi-channel pipettes of fixed and variable volumes, manual and automatic dispensers, burettes, dilutors and volumetric glassware.

The Suite of Quanta Software:

- Meets ISO 8655, DIN 12650, ISO 9000, ISO 17025 and GLP requirements
- Can be used with your current calibration protocols
- Connects to any balance supporting RS232 interface; controls and collects calibration data automatically from the balance
- Automatically converts mass into volume by using temperature and barometric pressure
- Provides automatic calculations on inaccuracy and imprecision
- Produces calibration reports and optional reports automatically
- Includes a database of 1500 pre-loaded test plans for all major manufacturers' pipettes and liquid handling instruments
- Provides calibration specifications in accordance with manufacturers' specifications, ISO 8655:2002 maximum permissible errors or user-defined specifications
- Maintains a complete service and calibration history, which can be easily accessed





Pipette calibration software
for the professional laboratory



Pipette service and calibration software
for the professional laboratory



Pipette service and calibration software
for the professional engineer

QuantaCal

Our entry level software. QuantaCal is designed primarily for end users that simply need to maintain a database of their pipettes, in order to check them against the manufactures' specifications, ISO 8655 or end user specifications, and to produce detailed calibration report that will satisfy general GLP | GMP requirements.

QuantaLab

Our midrange software. QuantaLab is a unique menu-based software system, incorporating QuantaCal as its base module and seven additional options

(see table). This software enables you to control the entire service and calibration process. It is ideally suited to smaller service companies or organisations that test and maintain their own pipettes in-house.

QuantaPro

The most comprehensive service and calibration software on the market. Developed for professional, high output pipette service companies or large multi-site end user accounts, this can be used to improve all aspects of the service process, from receipt to dispatch.

Quanta Software Modules	QuantaCal	QuantaLab	QuantaPro
Add Edit Pipettes	•	•	•
Personalised Information	•	•	•
RS232 Balance Support	•	•	•
Calibrate Pipettes	•	•	•
Add Edit Test Plans	•	•	•
Multiple Departments	•	•	•
Calibration Reports	•	•	•
Multiple Customer Entry	Optional	•	•
Booking In Dispatch System	Optional	•	•
Pipette Collection Requests ¹	Optional	•	•
Service Module	Optional	•	•
Spare Part Management ²	Optional	•	•
Due Overdue Pipettes Report	Optional	•	•
Service Calibration Recall	Optional	•	•
Email System	Optional	Optional	•
Quality Control	Optional	Optional	•
Report Editor	Optional	Optional	•
Management Reporting	Optional	Optional	•
Demo Equipment Management	Optional	Optional	•
Quote Generator	Optional	Optional	Optional
Mettler Multichannel Support	Optional	Optional	Optional
Accredited Certificates (ISO17025)	Optional	Optional	Optional
Multi-language Output	Optional	Optional	Optional

1) Pipette Collection Request is a sub-module of the Booking | Dispatch system. The Booking | Dispatch system must be in place before installing Pipette Collection Request.

2) Spare Parts Management is a sub-module of the Service Module. The Service Module must be in place before installing Spare Parts Management

Evaluation Version

Please download evaluation version from www.quantapro.net

Recommended System Requirements

- 2GHz Pentium P4 Processor
- Windows: 2000 | XP | Vista | Vista 64
- 256Mb RAM
- 100Mb Free Disk Space

How to Order?

- Please contact your local Sartorius representative for ordering information or
- Log on to www.quantapro.net and complete a contact form which will be forwarded to your nearest Sartorius representative

Online Support and Additional Information

Please visit www.quantapro.net for technical support and further information

Quanta software, developed by a professional pipette and service company, is truly the most comprehensive pipette service and calibration package available today.

 Ordering Information**Quanta**

Cat. No.	Item
760015	QuantaCal
760014	QuantaLab
760013	QuantaPro



■ Pipette Decontamination Procedure

Mechanical Pipettes (mLINE and Proline Plus)



1. Unscrew the tip ejector collar counter clockwise and remove it.



2. Unscrew the tip cone holder counter clockwise and carefully remove it along with the tip cone. Remove the Safe-Cone Filter if fitted.



3. Unscrew the piston counter clockwise from the pipette.



- Place the tip ejector collar, tip cone holder, tip cone and tip cone cylinder into a beaker containing 70% ethanol and leave for at least 30 minutes.
- After performing the procedure described above, remove the components from the beaker and rinse them with distilled water, then dry preferably with warm air, for at least an hour.
- Re-grease the piston as described in the instruction manual. Replace all components including new filter if fitted.

Electronic Pipettes (Picus and eLINE)



1. Unscrew the tip ejector collar counter clockwise and remove it.



2. Unscrew the tip cone holder counter clockwise and carefully remove the tip cone holder, tip cone and spring. Remove the Safe-Cone Filter if fitted.



3. Unscrew the piston counter clockwise from the pipette.



- Place the tip cone, tip cone holder, tip ejector collar, piston and spring into a beaker containing 70% ethanol and leave for at least 30 minutes.
- After performing the procedure described above, remove the components from the beaker and rinse them with distilled water, then dry preferably in warm air, for at least an hour.
- Re-grease the piston as described in the instruction manual. Replace all components including the new filter if fitted.

Note: When performing the decontamination procedure, as a matter of routine the o-ring seal should be checked for possible wear on every sixth occasion, and replaced if necessary. Calibration should also be checked at the same time.

Autoclaving Instructions



mLINE® and Proline Plus Mechanical Pipettes

The entire mLINE and Proline Plus mechanical pipette can be steam sterilized by autoclaving at 121°C (252°F), 1 bar (15 p.s.i.) for 20 minutes. The dispensing head of the multichannel pipettes must be unscrewed 360° counter clockwise before autoclaving.

- Remove the Safe-Cone Filter (if fitted)
- Put the pipette into the sterilisation bag and place it into the autoclave
- After autoclaving the pipette must be cooled down and left to dry overnight before use

It is recommended that you check the performance of the pipette after every autoclaving, and grease the piston | seal of the pipette after every 10th autoclaving.

Picus and eLINE® Electronic Pipettes Lower Parts

The dispensing head (tip ejector collar, tip cone holder, tip cone, spring and piston) of the single-channel and multi-channel models (except for multichannel 1200 µl) can be autoclaved (121°C, 1 bar for 20 minutes). These parts can be autoclaved as one unit or separately as individual parts. It is also possible to clean the parts and grease the piston prior to autoclaving.

- Remove the Safe-Cone Filter (if fitted).
- Put the dispensing head into the sterilisation bag and place it into the autoclave
- After autoclaving the parts must be cooled down and left to dry before use

It is recommended that you check the performance of the pipette after every autoclaving, and grease the piston | seal of the pipette after every 10th autoclaving.

Tips and Tip Boxes

- Place the bulk tips into the sterilisation bag and the tip tray as such in the autoclave
- Autoclave for 20 minutes at a temperature of 121°C under 1 bar (15 p.s.i.)
- Cool before use

Note:

- Excessive heat or length of time may damage the products. Never place the handle part of the Picus or eLINE into the autoclave
- The lower ends of multichannel pipettes are not interchangeable between 8 and 12-channel pipettes
- The cover of the tip tray should be closed during autoclaving

■ Troubleshooting Guide

Problem	Cause	Solution
Droplets left inside the tip	Unsuitable tip	Use original Sartorius tips
	Non-uniform wetting of the tip plastic	Attach new tip and pre-wet it
	Optiload not fully utilised and thus tip does not fit very well (All mLINE models and Proline Plus MC-models and 2ml, 5ml and 10ml SC-models have Optiload)	Pick up the tip so that it strikes the lower edge of the tip ejector collar
Leakage or pipetted volume too small	Tip incorrectly attached	Attach firmly
	Unsuitable tip	Use original Sartorius tips
	Tip is leaking and or Optiload not fully utilized	Replace a new tip or see above for Optiload
	Plunger movement not uniform, constant or balanced	It is very important that the plunger movement is slow and that this is always done in the same way during pipetting cycles. Volumes 1ml and greater: if plunger is released too quickly upon aspirating, it will affect the amount of liquid drawn.
Test results are incorrect and/or results are non-linear	Improper maintenance of lower parts	Clean and regrease lower parts as per instructions in this manual
	Piston or tip cone damaged/non-linear	Replace part with new one
	Dirt on stop surfaces/inside the machinery	Clean the stop surfaces/interior of the machinery as per instructions
	Uneven piston movement inconsistent pipetting technique/rhythm	It is very important that the plunger movement is slow and that this is always done in the same way during pipetting cycles. Volumes 1ml and greater: if the plunger is released too quickly upon aspirating, it will affect the amount of liquid drawn.
Display not sitting properly (in particular the lower edge of the display does not fit)	Calibration wheel (726066) not properly in place	Remove display and push calibration wheel down
Counter reading does not make sense and or counter feels very slack	The wings on the lead screw have come out of the grooves on the machinery	Replace the wings to the grooves
	The wings on the lead screw have been broken	Replace machinery with new equipment
	There is something broken inside the machinery	Replace machinery with new equipment
Multichannel pipettes: tip cone(s) does not draw liquid	Tip cone (bajonet connection) has accidentally come off	Refit tip cone
	Piston(s) is not connected to the piston support plate	Refit piston(s) so that they connect to the piston support plate properly
	Piston/tip cone damaged	Replace parts with new ones
	Improper maintenance	Maintain lower parts as per instructions

Problem	Cause	Solution
Pipette does not draw liquid at all	Magnet holder/magnet is disconnected from the lower part	Remove the dispensing head and refit it so that the magnet holder/magnet is properly connected to the dispensing head
Multichannel pipettes: Tip ejection does not work/ tip does not fit	Tip ejection bar has come off (the snap fit is disengaged)	Make sure that the tip ejection bar is properly connected
Electronic pipettes: Error in the display and motor is unable to start	Discharged battery/defective battery	Recharge battery/replace battery
	Actuator rod jammed	Clean and lubricate actuator rod
	Penetration of solvent vapours and thus actuator rod/tip ejection mechanism jammed	Clean tip ejection mechanism and clean/lubricate actuator rod
	Failure in the handle parts (upper body parts)	Check error messages
Electronic pipette can start but is unable to complete self test when switched on (error blinking on the display)	Discharged battery/defective battery	Recharge battery/replace battery
	Tip ejector mechanism jammed/contaminated	Clean tip ejection mechanism and clean/lubricate actuator rod
	Failure in the handle parts (upper body parts)	Check error messages
	Internal error has occurred	Check error messages and proceed accordingly
Tip ejector jammed or moves erratically	Ejector mechanism contaminated	Clean lower parts of pipette
	Ejector mechanism damaged	Replace damaged parts
Tip ejector feels slack	Ejector mechanism damaged	Replace damaged parts
Fading display and/or segment missing	Display damaged or incorrectly attached to the PCU-board/defective PCU board	Open handle cover and check the display
Autotest (a test programmed by service engineer) failed	Discharged battery/defective battery	Replace a battery
	Failure in the handle parts (upper body parts)	Check error counters and autotest counters and proceed accordingly
Pipette not charging	Pipette incorrectly positioned in the charging stand/carousel	Check the position of the pipette
	AC-adaptor damaged	Replace AC-adaptor
	Charging carousel/stand damaged	Open the bottom, check wires and replace the unit if needed
	Battery contacts inside the battery compartment damaged/oxidated/flattened	Clean battery contacts
	Battery damaged	Replace battery
Reduced operating time with fully charged batteries	Batteries damaged	Replace battery

■ Sales and Service Contacts

For Further Contacts, Visit www.sartorius.com

□ Europe

Germany

Sartorius Weighing Technology GmbH
Weender Landstrasse 94–108
37075 Goettingen

Phone +49.551.308.0
Fax +49.551.308.3289

France & Suisse Romande

Sartorius France S.A.S.
ZI Les Paluds
Avenue de Jouques – CS 31090
13781 Aubagne Cedex

Phone +33.442.845600
Fax +33.442.846545

Biohit France SAS

2, rue Antoine Laurent de Lavoisier
ZA de la Gaudrée
91410 Dourdan

Phone +33.170.625000
Fax +33 .1.64.59.76.39

Austria

Sartorius Austria GmbH
Franzosengraben 12
1030 Vienna

Phone +43.1.7965760.0
Fax +43.1.7965760.24

Belgium

Sartorius Belgium N.V.
Leuvensesteenweg, 248/B
1800 Vilvoorde

Phone +32.2.756.06.71
Fax +32.2.253.45.95

Hungary

Sartorius Hungária Kft.
Kagyló u. 5.
2092 Budakeszi

Phone +3623.457.227
Fax +3623.457.147

Ireland

Sartorius Ireland Ltd.
Unit 41, The Business Centre
Stadium Business Park
Ballycoolin Road
Dublin 11

Phone +353.1.8089050
Fax +353.1.8089388

Italy

Sartorius Italy S.r.l.
Via dell'Antella, 76/A
50012 Antella-Bagno a Ripoli (FI)

Phone +39.055.63.40.41
Fax +39.055.63.40.526

Poland

Sartorius Poland sp.z o.o.
ul. Wrzesinska 70
62-025 Kostrzyn

Phone +48.61.6473830
Fax +48.61.6473839

Russian Federation

LLC "Sartorius ICR"
and LLC "Sartorius Biohit"
Uralskaya str. 4, Lit. B
199155, Saint-Petersburg

Phone +7.812.327.5.327
Fax +7.812.327.5.323

Scandinavia, Finland & Baltics

Sartorius Biohit Liquid Handling Oy
Laippatie 1
00880 Helsinki

Phone +358.9.755.951
Fax +358.9.755.95.292

Spain & Portugal

Sartorius Spain S.A.
Offices in Madrid:
C/ Isabel Colbrand, 10–12, of. 70
28050 Madrid

Phone Spain +34.902.123.367
Phone Portugal +351.800.855.800
Fax Spain +34.91.358.84.85
Fax Portugal +351.800.855.799

Switzerland

Sartorius Mechatronics Switzerland AG
Ringstrasse 24a
8317 Tagelswangen (ZH)

Phone +41.44.746.50.00
Fax +41.44.746.50.50

U.K.

Sartorius UK Ltd.
Longmead Business Centre
Blenheim Road, Epsom
Surrey KT19 9QQ

Phone +44.1372.737159
Fax +44.1372.729927

□ America

USA

Sartorius Corporation
5 Orville Drive, Suite 200
Bohemia, NY 11716

Phone +1.631.254.4249
Toll-free +1.800.635.2906
Fax +1.631.254.4253

Argentina

Sartorius Argentina S.A.
Int. A. Ávalos 4251
B1605ECS Munro
Buenos Aires

Phone +54.11.4721.0505
Fax +54.11.4762.2333

Brazil

Sartorius do Brasil Ltda
Av. Dom Pedro I, 241
Bairro Vila Pires
Santo André
São Paulo
Cep 09110-001

Phone +55.11.4451.6226
Fax +55.11.4451.4369

Canada

Sartorius Canada Inc.
2179 Dunwin Drive #4
Mississauga, ON L5L 1X2

Phone +1.905.569.7977
Toll-Free +1.800.668.4234
Fax +1.905.569.7021

Mexico

Sartorius de México S.A. de C.V.
Circuito Circunvalación Poniente
No. 149
Ciudad Satélite
53100, Estado de México
México

Phone +52.5555.62.1102
Fax +52.5555.62.2942

□ Asia | Pacific

Australia

Sartorius Australia Pty. Ltd.
Unit 5, 7-11 Rodeo Drive
Dandenong South Vic 3175

Phone +61.3.8762.1800
Fax +61.3.8762.1828

China

Sartorius Scientific
Instruments (Beijing) Co., Ltd.
33 Yu An Road, Airport Industrial Park
Zone B,
Shunyi District, Beijing 101300,
P.R.China

Phone +86.10.8042.6300
Fax +86.10.8042.6486

Hong Kong

Sartorius Hong Kong Ltd.
Unit 1012, Lu Plaza
2 Wing Yip Street
Kwung Tong
Kowloon, Hong Kong

Phone +852.2774.2678
Fax +852.2766.3526

India

Sartorius Weighing India Pvt. Ltd.
#69/2-69/3, NH 48, Jakkasandra,
Nelamangala Tq
562 123 Bangalore, India

Phone +91.80.4350.5250
Fax +91.80.4350.5253

Japan

Sartorius Japan K.K.
4F Daiwa Shinagawa North Bldg.
8-11, Kita-Shinagawa 1-chome
Shinagawa-ku
Tokyo 140-0001

Phone +81.3.3740.5408
Fax +81.3.3740.5406

Malaysia

Sartorius Malaysia Sdn. Bhd
Lot L3-E-3B, Enterprise 4
Technology Park Malaysia
Bukit Jalil
57000 Kuala Lumpur, Malaysia

Phone +60.3.8996.0622
Fax +60.3.8996.0755

Singapore

Sartorius Singapore Pte. Ltd
1 Science Park Road,
The Capricorn, #05-08A,
Singapore Science Park II
Singapore 117528

Phone +65.6872.3966
Fax +65.6778.2494

South Korea

Sartorius Korea Ltd.
8th Floor, Solid Space B/D,
PanGyoYeok-Ro 220, BunDang-Gu
SeongNam-Si, GyeongGi-Do, 463-400

Phone +82.31.622.5700
Fax +82.31.622.5799

Thailand

Sartorius (Thailand) Co. Ltd.
129 Rama 9 Road,
Huaykwang
Bangkok 10310

Phone +66.2643.8361-6
Fax +66.2643.8367



Fapa s.a.s.
Via Pascoletto, 20
24040 Lallio (BG)
Tel. 035.6221219
Fax. 035.4372675
fapa@fapa.bg.it
www.fapa.bg.it

